


Practical Guide to Dressing Selection

For on-going chronic wound management make sure to undertake a full holistic assessment to address underlying conditions. Always work in a clean environment and aim to involve the patient and their caregivers.

Remove and check dressing



Does the dressing match the shape and depth of the wound?


- Shape** – A dressing with vertical absorption will not spread the wound fluid beyond the size of the wound, to help protect the surrounding skin.
- Depth** – Aim to prevent exudate pooling by selecting a dressing that maintains close contact with the wound bed.

Wound bed preparation



Healthy exudate

- Signs** – Healthy exudate (also called serous) is typically clear of color, thin or watery in consistency, and present in low to moderate amounts.
- Treatment** – At each dressing change, use saline or a wound cleanser such as Sea-Clens® to irrigate the wound, and gently pat dry before applying a dressing.



Suspected biofilm or infection

- Signs** – Thicker exudate (also called purulent) that is milky, yellow, green or brown, can be a sign of infection, especially in higher amounts and with a strong odor.
- Treatment** – In addition to irrigation, gentle physical debridement of the wound bed is recommended by using gauze or a pad with gentle circular motion in the wound bed.



Necrotic tissue such as slough or eschar

- Signs** – Dead or necrotic tissue (also called devitalized) often looks like thick slough or hard eschar (similar to a scab), and ranges from grey/white to yellow slough to black eschar.
- Treatment** – Refer to the dressing options below to help soften and absorb the necrotic tissue, but persistent slough or eschar may have to be referred for sharp or surgical debridement.
- Important** – Do not attempt to debride hard black eschar on the foot (below the ankle) or on diabetic patients without consulting a wound specialist.

Wound measurement



- Dimensions** – Measure the length and width of the wound bed and compare to previous measurements to see if the wound is healing.
- Depth** – Measure the depth of the wound at the deepest point and make a mental note to manage any depth when selecting a dressing.
- Undermining** – Check to see if there is any exposed wound bed underneath the skin by using a sterile probe, and document by imagining a clock on the wound bed and noting the hour to indicate where the wound is undermined (e.g. undermining at 8 o'clock).

Periwound skin



Macerated skin

- Signs** – When the surrounding skin (also called periwound) has been exposed to exudate, it could become macerated, looking soggy or soft and appearing white in color.
- Treatment** – Maceration can be prevented by using dressings that absorb wound fluid vertically, without spreading across. For skin that is already macerated, use a protective barrier film before applying the dressing.



Dry or flaky skin

- Signs** – Dry skin can be red, flaky or itchy (also called pruritus), sometimes with fine lines or cracks.
- Treatment** – Rehydrate the skin by using a moisture barrier cream, but make sure the cream has been completely absorbed before applying a dressing.

Dressing Selection

		Primary dressing	Secondary dressing
	Suspected biofilm/infection Treatment – Clean the wound bed with irrigation and gentle physical debridement. Treat with antimicrobial silver dressing. Systemic infections often cause fever-like symptoms and should be referred to a specialist immediately.	 Biatain® Alginate Ag Absorbs exudate and bacteria and contains antimicrobial silver	 Biatain® Silicone Foam Dressing Keeps Biatain Alginate Ag in place and offers additional absorption
	Slough or eschar Treatment – Thick slough or black eschar could be softened by adding moisture (e.g. hydrogel + secondary dressing). Persistent slough or eschar may need referral for sharp debridement.	 Triad® Hydrophilic Wound Dressing Alternative to hydrogel or collagenase ointment. Creates a moist environment to help soften slough or eschar and absorbs wound exudate and broken-down tissue	No secondary dressing is required, but gauze could be used to prevent brushing off
	Dry wounds (any depth) Treatment – For dry wounds with no wound fluid the primary goal is to add moisture to the wound bed to help the wound healing process.	 Woun'Dres® Collagen Hydrogel Rehydrates dry wounds and create a moist environment	 Biatain Silicone Foam Dressing Keeps Woun'Dres Collagen Hydrogel in place and helps to distribute it evenly
	Tunnelling or undermining Treatment – Dress any exposed wound bed underneath the skin with a primary filler dressing to fill any gaps and absorb excess exudate.	 Biatain® Alginate Primary filler dressing to manage cavities and undermining	 Biatain Silicone Foam Dressing Keeps Biatain Alginate in place and offers additional absorption
	Wounds with depth Treatment – Prevent exudate pooling in the wound bed by using a primary filler dressing to manage wounds with depth.	 Biatain Alginate Creates close contact with the wound bed to prevent pooling	 Biatain Silicone Foam Dressing Keeps Biatain Alginate in place and offers additional absorption
	Superficial wounds Treatment – Partial-thickness wounds or wound with minimal depth can often be managed by a silicone foam dressing.	 Biatain Silicone Foam Dressing Conforms to the wound bed to prevent any small gaps, and absorbs wound fluid vertically without spreading across onto surrounding skin	No secondary dressing is required, Biatain Silicone Foam Dressing acts as a primary and secondary dressing in one
	Difficult-to-dress Treatment – Aim to keep all wounds covered to create a moist wound healing environment, even in difficult areas such as gluteal cleft, wounds near incontinence, wounds with macerated skin.	 Triad Hydrophilic Wound Dressing Sterile coating for wet wounds or wet skin, or for difficult-to-dress areas. Keeps wounds covered and absorbs medium levels of wound fluid	No secondary dressing is required, but gauze could be used to prevent brushing off
	Prevention or stage 1 Treatment – Foam dressing to help manage pressure, shear/friction, and moisture on the skin – in combination with your local prevention protocol.	 Biatain Silicone Foam Dressing Helps to address risk factors by redistributing pressure, reducing friction and shear, and managing microclimate	No secondary dressing is required