

Bladder & Bowel Management Options

SpeediCath® intermittent catheters available in smaller French sizes & Peristeen® Plus Transanal Irrigation System



All SpeediCath catheters feature **Triple Action Coating Technology** to protect the urethra and reduce the risk of urinary tract infections.^{1*} There are SpeediCath options available in a variety of lengths and French sizes to meet your patients' specific needs.

Available in smaller French sizes



SpeediCath® Standard Boy

10" length

Code	Size
28606	● 6 FR
28608	● 8 FR
28610	● 10 FR
28612	○ 12 FR



SpeediCath® Standard

14" length

Code	Size
28408	● 8 FR
28410	● 10 FR
28412	○ 12 FR



SpeediCath® Soft

13" length

Code	Size
27010	● 10 FR
27012	○ 12 FR
27014	● 14 FR



SpeediCath® Flex Coudé Pro

13" length

Code	Size
20021	● 10 FR - Pocket
20022	○ 12 FR - Pocket
20024	● 14 FR - Pocket
20010	● 10 FR
20012	○ 12 FR
20014	● 14 FR



SpeediCath® Compact Set

13.2" length

Code	Size
28422	○ 12/18 FR



SpeediCath® Flex Set

13" length

Code	Size
28931	● 10 FR
28932	○ 12 FR
28934	● 14 FR



*Compared to uncoated catheters.

1. Rognoni C, Tarricone R. Intermittent catheterisation with hydrophilic and non-hydrophilic urinary catheters: systematic review. *BMJ* 2007;335:e1258. PMID: PMC5225586.

on Coating Technology



Stays bonded, stays smooth, and stays hydrated for the protection of the urethra*



SpeediCath® Standard

6" length

Code	Size
28506	● 6 FR
28508	● 8 FR
28510	● 10 FR
28512	○ 12 FR



SpeediCath® Compact Set

3.5" length

Code	Size
28520	● 10 FR
28522	○ 12 FR
28524	● 14 FR



SpeediCath® Standard Pediatric

6" length

Code	Size
28706	● 6 FR
28708	● 8 FR
28710	● 10 FR



SpeediCath® Compact Female and Compact Plus

2.75" length

Code	Size
28578	● 8 FR
28580	● 10 FR
28582	○ 12 FR
28584	● 14 FR



3.5" length

Code	Size
28810	● 10 FR
28812	○ 12 FR
28814	● 14 FR



Peristeen® Plus

Peristeen® Plus Transanal Irrigation System

For predictable bowel
management



Transanal irrigation (TAI) is an option for children and adolescents with neurogenic bowel dysfunction that can help them:

- Take control - decide when and where to empty the bowels
- Help to prevent constipation and fecal incontinence
- Minimize the likelihood of involuntary release of stool for up to two days
- Improve their quality of life²

Peristeen® Plus Transanal Irrigation

Code

29152	Peristeen Plus system without catheters	1 control unit / 1 water bag (with lid) / 2 straps / 1 tube
29142	Peristeen Plus Accessory unit Regular	15 rectal balloon catheters* Regular / 1 water bag
29149	Peristeen Plus Accessory unit Small	15 rectal balloon catheters* Small / 1 water bag

Peristeen Plus System

Peristeen Plus is a prescription-only medical device intended to instill water into the colon through a rectal catheter – which incorporates an inflatable balloon – inserted into the rectum to promote evacuation of the contents of the lower colon.

2. Christensen P, et al. A randomized, controlled trial of transanal irrigation versus conservative bowel management in spinal cord-injured patients. Supported by Coloplast. Gastroenterology 2006;131:738-747, n=87.



Coloplast® Care is a product and lifestyle support program for intermittent catheter and Peristeen® Plus Transanal Irrigation System users



Coloplast Care Advisor
Individualized follow-up by phone



Product Samples*
and Tools




Clinician-validated Product and
Lifestyle Education

Coloplast Care guides your patients from
confusion to *confidence*

Get started today!

Enroll patients in Coloplast Care, select samples,* and send their order to a supplier that works with their insurance – all within one platform.



To learn more or request **Coloplast Care Connect** with **Epic®** integration† for your practice, please contact your Coloplast representative or visit www.coloplastcareconnect.us

Scan the QR code with your smartphone camera and select “Start Here” to enroll your first patient



Information from Coloplast® Care is for educational purposes only. It is not intended to substitute for professional medical advice and should not be interpreted to contain treatment recommendations. You should rely on the healthcare professional who knows your individual history for personal medical advice and diagnosis.

*Limitations apply. Intermittent catheter samples available with valid prescription only.

Epic® is a registered trademark of Epic Systems Corporation.

†Epic integration is optional and is not required to use the Coloplast Care Connect portal.

Prior to use, refer to product labeling for complete product instructions for use, contraindications, warnings and precautions.

SpeediCath® catheters are prescribed for use by patients who require bladder drainage due to urine retention or post void residual volume (PVR). Before use, carefully read all of the instructions. Call your doctor if you think you have a UTI or can't pass the catheter into the bladder. For more information regarding risks, potential complications and product support, call Coloplast Corp. at 1-866-226-6362 and/ or consult the company website at www.coloplast.us.

Indication: The Peristeen® Plus Transanal Irrigation is intended to instill water into the colon through a rectal catheter, which incorporates an inflatable balloon, inserted into the rectum to promote evacuation of the contents of the lower colon. The Peristeen Plus Transanal Irrigation System is indicated for use by children (2 years - <12 years old), adolescent (12 years - < 18 years old), transitional adolescent (18 - <21 years old) and adult patients with neurogenic bowel dysfunction who suffer from fecal incontinence, chronic constipation, and/or time-consuming bowel management procedures.

Contraindications: Peristeen Plus Transanal Irrigation must not be used in the following situations: known anal or colorectal stenosis, colorectal cancer (active/recurrent), radiotherapy to the pelvis, and recent abdomino-perineal surgery, active inflammatory bowel disease, acute diverticulitis, severe diverticulitis, previous diverticulitis and diverticular abscess, chronic symptomatic diverticular disease, within 3 months of abdominal, anal or colorectal surgery, within 4 weeks of endoscopic polypectomy, recent colonic biopsy, recent endoscopic mucosal resection (EMR) and recent endoscopic sub-mucosal dissection (ESD), severe autonomic dysreflexia, ischemic colitis, during Spinal Cord Shock Phase, complex diverticular disease, in patients who are pregnant and have not used the system before.** Since the list is not exhaustive, the healthcare professional should always consider individual patient factors as well.

**If the patient is pregnant and has never used transanal irrigation before, they should not start the irrigation procedure during pregnancy.

Warning: Peristeen Plus Transanal Irrigation procedure should always be carried out with caution. Bowel perforation is an extremely rare, but serious and potentially lethal complication to transanal irrigation and will require immediate admission to a hospital, often requiring surgery. Seek immediate medical assistance if during or after Peristeen Plus Transanal Irrigation procedure you experience any of the following: severe or sustained abdominal pain or back pain, especially if combined with fever; severe or sustained anal bleeding; blood in the feces; rectal pain, and abdominal pain. See the device manual for complete user instructions, contraindications, warnings, cautions. For further information, call Coloplast Corp at 1-855-605-7594 or consult the company website at www.coloplast.us. Caution: Federal law restricts this device to sale by or on the order of a physician.

Information from Coloplast® Care is for educational purposes only. It is not intended to substitute for professional medical advice and should not be interpreted to contain treatment recommendations. You should rely on the healthcare professional who knows your individual history for personal medical advice and diagnosis.