

SpeediCath[®] Compact Set

Discreet and easy to use

All-in-one
catheter



SpeediCath[®] Compact Set

An all-in-one catheter and bag solution



reddot design award
winner 2013

SpeediCath® Compact Set is part of the award winning SpeediCath range.



Available in
male and
female lengths



Male



Female

All-in-one

Integrated urine bag

No need to connect a urine bag, get to a toilet or store a bulky product.

Instantly ready-to-use

No need to add water before use, thanks to its smooth hydrophilic coating that also minimizes friction.

Simple

Easy opening and re-closing

Can be opened and closed in a single twist.

Easy grip handle

Hold without touching the pre-lubricated part of the catheter.

Discreet

Designed for discretion

Easy to bring along and store with its discreet hard case.

Product details

Female

Code	Size	Units	Description	HCPCS*
28520	● 10 FR	20	3.5" length	A4353
28522	○ 12 FR	20	3.5" length	A4353
28524	● 14 FR	20	3.5" length	A4353

Male

28422	○ 12/18 FR	20	13.2" length includes funnel	A4353
-------	------------	----	------------------------------	-------

* Reimbursement Disclaimer: Coloplast Corp. provides this information for your general reference and related to the reimbursement of Coloplast products only. Reimbursement, coverage and payment policies can vary from one insurer and region to another, and may change over time. Coloplast does not guarantee coverage or payment of products.

SpeediCath® Compact Set catheters are prescribed for use by male or females who require bladder drainage due to urine retention or post void residual volume (PVR). Before use, carefully read all of the instructions. Call your doctor if you think you have a UTI or can't pass the catheter into the bladder. For more information regarding risks, potential complications and product support, call Coloplast Corp. at 1-866-226-6362 and/or consult the company website at www.coloplast.us