Ostomy
Product Guide



Mike
User I SenSura® Mio Convex Urostomy Pouch

Welcome to Coloplast

Supporting Products

Irrigation

2-Piece Click

2-Piece Flex

1-Piece

Fistula/Wound Management

Pediatric

Services

Coloplast

Effective November 2026

Welcome to Coloplast

Coloplast was founded on passion, ambition and commitment. We were born from a nurse's wish to help her sister and the skills of an engineer. Guided by empathy, our mission is to make life easier for people with intimate healthcare needs.

Over decades, we have helped millions of people live more independent lives and continue to do so through innovative products and services.

The right fit is key to reducing leakage

Leakage affects the majority of ostomy patients[1]



76%
of ostomates experience leakage[1]

66%
of nurses consider leakage to be a main reason for peristomal skin complications[2]

BodyFit Technology is an elastic adhesive that is unique to all **SenSura® Mio** barriers

- **Adheres securely** to the body to protect the skin
- **Flexible** to maintain a secure fit during body movement
- **Flows** into creases and folds to prevent leakage

BodyFit Technology®



The above image is the SenSura Mio barrier conforming to the curves of the hand and flexing with movement

1. Claessens I, Probert R, et. al. The Ostomy Life Study: the everyday challenges faced by people living with a stoma in a snapshot. *Gastrointestinal Nurs.* 2015;13(5):18-25.
2. Fellows J, Voegeli D, Håkan-Bloch J, et. al. Multinational survey on living with an ostomy: prevalence and impact of peristomal skin complications. *British Journal of Nursing.* 2021;30(16):947-S44.

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Because every *fit* matters

Finding the right fit is essential to help reduce leakage and skin complications



Flat

Features a unique elastic adhesive that follows body shape and movement.



Convex

Provides stability without compromising flexibility. Available in Soft, Light, and Deep to fit a wide range of patient needs.





Convex Flip

The only barrier specifically designed for outward curves. Just flip to fit!

Peristomal body profiles vary greatly


A multi-site US acute care study showed SenSura® Mio pouching systems had 41% lower risk of leakage compared to users of Hollister® CeraPlus™ pouching systems*[4]

SenSura Mio offers a wide range of options to achieve the right fit to individual bodies and prevent the cycle of leakage and skin complications. What's more, our convex and convex flip products have each been proven to reduce leakage.[1-3]

Features of SenSura Mio Convex barriers

Designed with double layer, elastic adhesive

- Maintains secure contact to the skin
- Provides a flexible fit during movement



Convex FlexShell™ with integrated flexlines

- Provides stability without compromising flexibility
- Special stability and flexibility zones maintain secure contact and peristomal support, even during movement
- Helps with stoma protrusion when needed

1. SenSura Mio Convex Soft compared to user's own flat and reference soft convex product (Pelican Select Convex); n=38. Kruse TM, Størling ZM. Considering the benefits of a new stoma appliance: A clinical trial. *Br Jour Nurs.* 2015;24(22):S12-18.
2. SenSura Mio Convex Light and Deep compared to SenSura Convex Light; n=129. Walker H, Hopkins G, Waller M, Størling ZM. Raising the bar: New flexible convex ostomy appliance — a randomised controlled trial. *World Council of Enterostomal Therapists Journal.* 2015;35(1):S6-11.
3. SenSura Mio Convex Flip (referred to as 'concave') compared to a 2-piece flat pouching system; n=52. Rolfsen T; Vestergaard M, Hansen MF, Boisen EB, Dambæk MR. Body fit with a pouching system with concave contour for people with an outward peristomal body profile: Effects on leakage, wear time, and quality of life: a randomized controlled cross-over trial. *J Wound Ostomy Continence Nurs.* 2024;51(4):303-311.
4. Based on a retrospective study conducted at 2 affiliated university hospitals; n=214 patient charts with surgery dates ranging from February 2022 to June 2023. Patient events were charted by the WOC nurse at each pouch change; review captured up to 5 pouch changes or 2 weeks of data. Statistical analysis controlled for key differences between sites and found predicted probability of leakage for SenSura Mio was 15.1% and Hollister® CeraPlus™* pouching system was 25.6%, corresponding to a 40.8% lower risk (p=0.011). Zamarripa, C. et al. Pouching system leakage and peristomal skin complications following ostomy surgery in the immediate postoperative period: A retrospective review. *J Wound Ostomy Continence Nurs.* 2024;51(6):478-483.

™* All trademarks are the property of their respective owners.

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

SenSura® Mio Convex

A full range of options – with proven leakage reduction – to meet the needs of all of your patients

SenSura® Mio Convex Soft - 6mm



High flexibility, high compressibility, gentle support

Patient profile:
- Distended or firm abdomen
- Small creases
- Protruded stoma with leakage
- Postoperative abdomen




SenSura® Mio Convex Light - 6mm



High flexibility, medium compressibility, more stable

Patient profile:
- Crease/fold
- Stoma opening not centered
- Inward shape
- Soft abdomen or tissue
- Postoperative abdomen



SenSura® Mio Convex Deep - 9mm



High flexibility, low compressibility (most stable). Deepest on the market

Patient profile:
- Large crease/fold
- Stoma is below skin level or retracted
- Stoma opening not centered
- Inward shape
- Postoperative abdomen



New!



SenSura® Mio Urostomy Pouch with *Twist Outlet*

Redesigned faucet-style outlet opens and closes with a simple twist

Enhanced usability and hygiene with a winged-grip design and hygienic emptying experience* that helps keep urine off patients' hands.

Secure, leak-tight twist-lock outlet designed to reduce urine drips and odor,* with clear confirmation of closure and minimal risk of accidental opening.

Smooth, controlled drainage enabled by a wide-lumen outlet for efficient emptying.*



Twist to open
The outlet is open when the tap is turned counterclockwise



Twist to close
The outlet is closed when aligned with the pouch

Empower patient independence in everyday life with the ease of use and hygienic experience* of SenSura Mio urostomy pouch with new twist outlet.



Preferred by users and Healthcare Professionals when compared to Hollister® and ConvaTec®*

* Simulated use, Coloplast Data-on-File (09/2023), n=30

All trademarks and copyrights are the property of their respective owners.


Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

A new adapter for urostomy night drainage bags

SenSura® Mio urostomy pouches come with a new adapter so they can connect to your choice of universal bedside drainage bag.**

Always use the adapter to connect to a bedside night drainage bag

With the adapter you can connect the SenSura Mio urostomy pouches outlet to any universal bedside drainage bag,** including those from Coloplast. Please note that you can no longer use bedside night drainage bags without the adapter included in the box.




**When using the urostomy adapter

SenSura® Mio Urostomy 2-Piece Click now with *Improved Coupling[2] and New Twist Outlet*

SenSura Mio Urostomy 2-Piece Click coupling with added security

- Gentle, secure and leak-tight†
- Preferred over other coupling systems†
- Application is easier to learn and teach[1,2] - **Place, Click, Lock**



1. Coloplast Data-on-File, Product Evaluation by HCPs (11/2025), n=345.
2. Compared to the original SenSura Mio Click.
†Simulated use, Coloplast Data-on-File (09/2023), n=30

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Ostomy care for skin, body and mind

Supporting your care of every patient
Our complete portfolio of 1-piece and 2-piece pouching systems can meet the unique needs of your patients.

2-piece Click pouching systems
Mechanical coupling with a locking mechanism for security you can hear



2-piece Flex pouching systems
Adhesive coupling provides the convenience of a 2-piece with the flexibility of a 1-piece



1-piece pouching systems
For a low profile and ultimate flexibility



Innovative pediatric solutions *for little fighters*


SenSura® Mio Baby
For babies 30 weeks' gestation to 6 months





- Specifically designed to **meet the needs of premature babies and infants** living with a stoma
- Includes a **refeeding window** for easy access to refeed or expel gas
- **Skin-friendly adhesive** provides long wear time for less disturbances

SenSura® Mio Kids
For kids 6 months to 4 years


Henry
SenSura® Mio Kids user



- **Star-shaped elastic adhesive** has special petals to hug kids' small, round bellies
- Designed with Bodyfit Technology to provide a **flexible fit** to allow children to move freely
- Options available in 1-piece and 2-piece pouching system to **step up as kids grow**

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Innovative supporting products ***to personalize fit***

Brava® Supporting Products

Our barrier and pouching systems in combination with Brava® Supporting Products provide a personalized pouching solution to help ***reduce leakage*** and ***help reduce skin issues.***

Protect the skin


Brava® Adhesive Remover Spray and Wipes


Brava® Protective Sheet


Brava® Skin Barrier Spray and Wipes


Brava® Powder


Brava® Skin Cleanser Wipes

Secure against leakage


Brava® Protective Seal


Brava® Protective Seal Convex


Brava® Strip Paste


Brava® Paste

Support with confidence


Brava® Elastic Barrier Strips


Brava® Lubricating Deodorant


Brava® Belt for SenSura® Mio


Brava® Ostomy Support Belt

For help eliminating airborne odors, check out Hex-On Odor Antagonist!




Lucie
User | SenSura® Mio and Brava®

Brava® Supporting Products

Personalized solutions for everyday life with an ostomy

Leakage and skin irritation are two of the most common problems faced by people with an ostomy, and they can carry severe consequences.

76% of ostomates experience leakage[1]

88% of ostomates experience skin issues[2]

95% of all people with an ostomy use supporting products[3]

3.7 supporting products are used on average[3]

1. Claessens I, Probert R, et. al. The Ostomy Life Study: the everyday challenges faced by people living with a stoma in a snapshot. *Gastrointestinal Nursing.* 2015;13(5):18-25.
2. Fellows J, et al. Multinational survey on living with an ostomy: prevalence and impact of peristomal skin complications. *Br J Nurs.* 2021;30(16):S22-S30.
3. Ostomy Life Study 2019, Data on file.

Brava® Supporting Products address three essential needs for everyday life with an ostomy

Protect the skin

Skin issues

Support healthy peristomal skin and help address common skin challenges.



Adhesive Remover | Skin Barrier | Protective Sheet | Powder | Skin Cleanser

Secure against leakage

Leakage

Create a secure seal to help prevent leakage and protect the skin around the stoma.



Protective Seal | Protective Seal Convex | Paste | Strip Paste

Support with confidence

Security

Provide added comfort, discretion, and confidence to help people with an ostomy feel supported throughout everyday life.



Elastic Barrier Strips | Lubricating Deodorant | Belt | Ostomy Support Belt



For help eliminating airborne odors check out Hex-On Odor Antagonist!

Brava® Supporting Products

Personalized solutions for everyday life with an ostomy

Brava® Elastic Barrier Strips

Skin-friendly strips are used to "picture frame" around an ostomy barrier for extra security.

- Elastic to follow body contours
- Supports longer wear-time by reducing lifting and rolling of the barrier edges
- Reduces unplanned changes
- Absorbs moisture and even sweat while being active
- Skin-friendly and secure alternative to tape
- Not made with natural rubber latex

Without stretching, apply strips after the barrier has been placed



Brava® Elastic Barrier Strips - Curved

For round or oval barrier shapes

Code	Size	Units	HCPCS*
120700**	5½" strips	20 strips	A4362

**2 strips = 1 billable unit of A4362



Brava® Elastic Barrier Strips - Y Shape

Wings provide additional support even in areas like hernia, bulges or deep folds

Code	Size	Units	HCPCS*
120721**	5¾" strips	30 strips	A5121

**3 strips = 1 billable unit of A5121



Brava® Elastic Barrier Strips - Straight

For square or rectangular barrier shapes

Code	Size	Units	HCPCS*
120740**	5¼" strips	20 strips	A4362

**2 strips = 1 billable unit of A4362



Brava® Elastic Barrier Strips - XL

60% larger adhesive than original Curved for increased security & support around the barrier

Code	Size	Units	HCPCS*
12076**	6¾" strips	30 strips	A5121

**3 strips = 1 billable unit of A5121


*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Brava® Protective Seal

Proven Leakage Protection

Brava® Protective Seal is **proven** to significantly reduce leakage by creating a **long-lasting** and **personalized** seal between the barrier and skin.



Brava® Protective Seal



Moldable barrier ring designed to adhere securely and protect the skin around the stoma.

- Easy to shape, apply, and removes cleanly in one piece
- Mold to fit any shape of stoma, including oval
- Soft, flexible material adapts to body shape as the wearer bends and stretches
- Unique polymer formulation is designed to hold up against even the most aggressive output
- Alcohol-free and sting-free
- Multiple sizes for a personalized fit
- Not made with natural rubber latex

Units 10, HCPCS* A4385

Stoma size	Code		Seal width
	Thin (2.5 mm)	Thick (4.2 mm)	
¾" (18 mm)	**12035****	**12045****	Standard (⅝", 15 mm)
	12036	**12046**	Wide (¾", 20 mm)
	12038	**12048**	Wide (⅞", 23 mm)
	12032	-	Extra wide (1⅛", 29 mm)
1 ⅛" (27 mm)	**12037****	**12047**	Standard (⅝", 15 mm)
	12033	-	Extra wide (1", 25 mm)
1 ⅜" (34 mm)	**12039**	**12049**	Standard (⅝", 15 mm)

**Denotes seals that fit most patients

Brava® Protective Seal Convex



Moldable barrier ring with convexity to provide support for flush or retracted stomas and patients with deep creases, folds or scars around their stoma.

- Sustains convexity after molding to fit individual stoma shapes, including oval and loop
- High-tack adhesive that flows quickly into skin, filling even the smallest creases
- Erosion resistant composition
- High flexibility to maintain seal during body movement
- Steep sides of convex seal are designed to provide convexity
- Thin edges fit snugly around peristomal area
- Not made with natural rubber latex

Units 10, HCPCS* A4411

Stoma size	Code	Thickness
¾" (20 mm)	**12090****	8.2 mm
1" (25 mm)	**12091**	8.2 mm
1 3/16" (30 mm)	**12092****	8.2 mm
1 ⅜" (35 mm)	**12093**	8.2 mm
1 9/16" (40 mm)	**12094****	8.2 mm
1 ¾" (45 mm)	**12095**	8.2 mm

**Denotes seals that fit most patients

Six must-have seals

Brava® Protective Seals are all designed to stretch, mold, and fill creases to even out the peristomal area. Below we've outlined **six seals** that meet most clinical needs.

Seals that fit most patients

For smaller stomas

Brava Protective Seal #12035
HCPCS* A4385

- Thin (2.5 mm) moldable ring with a ¾" (18 mm) opening
- Stretches to fit the most common sizes of stomas



For larger stomas

Brava Protective Seal #12037
HCPCS* A4385

- Thin (2.5 mm) moldable ring with a 1 ⅛" (27 mm) opening
- Stretches to fit larger stomas





For challenging peristomal areas

Brava Protective Seal #12045
HCPCS* A4385

- Thick (4.2 mm) moldable ring with a ¾" (18 mm) opening
- Fills creases, folds, and scars next to the stoma
- Can shape and use like strip paste to caulk and even out peristomal area prior to sealing area between stoma and barrier with additional ring fit to stoma size**





***May use more than 1 seal on a single patient**

Seals for patients that need convexity

For flush or retracted stomas or deep creases, folds or scars

Brava® Protective Seal Convex
HCPCS* A4411


#12090: ¾" (20 mm) opening
#12092: 1 3/16" (30 mm) opening
#12094: 1 9/16" (40 mm) opening

- Select convex seal closest to stoma size, then stretch to fit
- Moldable convex ring with 8.2 mm of convexity
- Convex shape protrudes stoma helping output go into pouch
- Sustains convexity even after being stretched and molded






A versatile seal for a personalized fit



Use as a ring to **support the stoma**



Shape and use like caulk to **fill folds, creases, and scars** to even out the peristomal area



Seal and fill gaps between the barrier and the stoma

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Brava® Moldable Ring



Legacy moldable barrier ring fills uneven skin surfaces to create a durable seal between the stoma and ostomy barrier to protect against leakage.

- Helps keep skin healthy by protecting it from output and absorbing moisture from the skin
- Alcohol-free and sting-free
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
120307	2 mm (thin)	10	A4385
120427	4.2 mm (thick)	10	A4385

Brava® Strip Paste



Easy to shape and mold paste strips, fill uneven surfaces for a flat pouching surface.

- Protects peristomal skin from stoma output
- Contains pectin
- Alcohol-free and sting-free
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
26555	2 oz. total package (0.2 oz. per strip)	10 strips	A4406

Brava® Paste



Paste is sting-free and is used to fill uneven skin surfaces for a flat pouching surface.

- Use like caulk to fill creases, folds and scars
- Creates a tight seal between peristomal skin and the ostomy/fistula appliance
- A sting-free alternative to conventional alcohol-based paste
- **No** alcohol = no need to "set"; **contains** pectin
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
12050	2 oz.	1	A4406

Brava® Skin Cleanser Wipes



The wipe that cleans and cares. A convenient way to cleanse and maintain healthy skin when changing the ostomy pouching system.

- Does not contain alcohol or perfumes
- No residue left on the skin that could affect barrier adhesion
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
12080	1 wipe	15 wipes/pack	A9286**

Supporting Products

Brava® Lubricating Deodorant



A pouch lubricant and deodorizer in one. Add inside pouch to neutralize odor and create a slick surface, allowing output to slide to the bottom of the pouch for easier emptying.

- Powder fresh smell
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
12060	.25 fl. oz. (7.5 mL) sachet	20 sachets	A4394
12061	8 oz. bottle	16 bottles	A4394

Brava® Adhesive Remover Spray



Quick and easy removal of adhesives.

- Designed to lift the edges of the barrier for ease in removal
- Dries within seconds without leaving an oily residue behind
- Alcohol-free, silicone-based, and **sting-free**
- Spray at any angle
- Spray can be completely emptied, reducing waste
- It does not change adhesion of next pouching system
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
120105	1.7 oz. (50 mL)	1	A4455**

Brava® Adhesive Remover Wipes



Wipe away adhesive residue.

- Skin-friendly wipes contain silicone to be gentle on the skin and allows for easy removal of the barrier
- Dries within seconds without leaving an oily residue behind
- Alcohol-free, silicone-based, and **sting-free**
- It does not change adhesion of next pouching system
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
120115	1 wipe	30 wipes/box	A4456**

Brava® Skin Barrier Spray



Creates an invisible skin friendly layer to protect the skin from output and adhesives.

- Dries matte and within seconds without building a thick or tacky layer
- Alcohol-free, silicone-based, and **sting-free**
- Improves skin condition by reducing skin problems related to output and adhesives
- Adhesion of next pouching system is unaffected
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
120205	1.7 oz. (50 mL)	1	A4369**

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

**Product has not been through PDAC verification

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Brava® Skin Barrier Wipes



Convenient wipe creates an invisible skin friendly layer to protect the skin from output and adhesives.

- Improves skin condition by reducing skin problems related to output and adhesives
- Adhesion of next pouching system is unaffected
- Dries matte and within seconds without building a thick or tacky layer
- Alcohol-free, silicone-based, and **sting-free**
- Not made with natural rubber latex

Tip: *Wring the wipe like a wet rag to expel the liquid, then dab to disperse on skin*

Code	Size	Units	HCPCS*
120215	1 wipe	30 wipes/box	A5120**

Brava® Protective Sheet



Skin friendly hydrocolloid skin barrier that protects the skin from stoma output.

- Adheres to moist skin
- Absorbs moisture and minimizes the risk of maceration
- Easy to cut to fit stoma or desired shape
- Provides a ready surface for good adhesion
- Not made with natural rubber latex

Without stretching, apply sheet before adding a seal or the pouching system

Code	Size	Units	HCPCS*
32105	4 x 4" (10 x 10 cm)	10	A4362
32155	6 x 6" (15 x 15 cm)	5	A5121
32205	8 x 8" (20 x 20 cm)	5	A5122

Brava® Powder



Non-sterile powder for moisture absorption to reduce risk of skin irritation.

- May use in conjunction with Brava® Skin Barrier Spray/Wipes or Prep™ Protective Skin Barrier to "crust", thereby protecting moist skin
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
19075	1 oz. (25 g)	16	A4371

Hex-On® Odor Antagonist



Air spray tackles and controls airborne odors.

- Non-aerosol spray
- Available in both Fresh Linen and Fragrance Free scent
- Not made with natural rubber latex

Code	Product	Size	Units
7583	Fresh Linen	2 fl. oz./59 mL	12
67589	Fragrance Free	2 fl. oz./59 mL	12

Ostomy Scissors



Curved ostomy scissors.

- Not made with natural rubber latex

Code	Size	Units	HCPCS*
95050	One size	30/bx	N/A

Supporting Products

Urostomy Adapter



New

- Urostomy adapter - connect the new SenSura® Mio urostomy pouch's outlet to any universal bedside drainage bag
- Easy clip-on mechanism - secures to the outlet with a gentle pinch enabling connection to bedside drainage bag

Code	Units	HCPCS*
18790*	5 adapters	NA

*Purchase available for hospital / acute settings only.
Patients who need additional adapters can look to request for samples by calling Coloplast Care 1-877-858-2656

Coloplast® Skin Barrier Rings



Soft and flexible legacy barrier rings provide skin protection for stomas, fistulas, and draining wounds.

- Not moldable
- Non-pectin based
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
2310	⅜" (10 mm)	30	A4404
2315	⅗" (15 mm)	30	A4404
2320	¾" (20 mm)	30	A4404
2325	1" (25 mm)	30	A4404
2330	1⅛" (30 mm)	30	A4404
2340	1⅗" (40 mm)	30	A4404
2350	2" (50 mm)	30	A4404

Prep™ Protective Skin Barrier Wipes



Clear film forms a breathable tacky layer protecting against skin irritants such as output and adhesives.

- No need to remove between applications
- Contains alcohol
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
2041	54 wipes/pack	12 packs/box	A5120

Filtrodor® Pouch Filters



Charcoal-activated, deodorizing filter with self-sealing foam center.

- Adheres directly to pouch
- In sheets of 10 filters with a needle to create perforation for airflow
- Box of 50 filters includes 5 needles
- Not made with natural rubber latex

Code	Units	HCPCS*
0509	50 *each / box contains 5 packs of 10*	A4368

Coloplast® Pouch Clamp



Simple closure clamp for drainable pouch, is easy to manipulate, even for those with low dexterity.

- Not made with natural rubber latex

Code	Units	HCPCS*
9500	20	A4363

*For important information regarding the use of HCPCS codes listed and code definitions see page 124
**Product has not been through PDAC verification

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Brava® Ostomy Support Belt



This belt relieves the heavy sensation of hernias and is designed to prevent hernias post-operatively.

- Comfortable and breathable fabric
- Stretchable to fit body contours
- Intuitive pocket closure
- Includes a silicone grip to prevent sliding and rolling
- Cut to fit panel can be used to have pouch on outside
- Not made with natural rubber latex
- All sizes are an 8" width

Tip: *Measure before ordering, sizes specific to belt*

Code	Size	Units	HCPCS*
12003	S (29-32")	1	A4396
12004	M (31-35")	1	A4396
12005	L (34-39")	1	A4396
12006	XL (37-42")	1	A4396
12007	2XL (40-46")	1	A4396
12008	3XL (44-51")	1	A4396
12009	4XL (50-58")	1	A4396

Brava® Belt for SenSura® Mio



Belt is designed with double hooks for SenSura® Mio barriers to provide extra stability across the entire peristomal area.

- Offers extra security for people with a flush or retracted stoma by supporting a convex barrier and allowing the stoma to protrude
- Adjustable, comfortable and washable
- Neutral gray color blends in under clothing
- Four tabs connect belt to SenSura Mio products with belt tabs
- Not made with natural rubber latex

Tip: *Measure before ordering*

Code	Size	Units	HCPCS*
4237	Up to 40" length	15	A4367
4247	Up to 61" length	15	A4367

Brava® Belt



Designed with single hooks to be compatible with most conventional ostomy products with belt tabs. The belt provides extra stability across the entire peristomal area.

- Offers extra security for people with a flush or retracted stoma by supporting a convex barrier and allowing the stoma to protrude
- Adjustable, comfortable and washable
- Neutral gray color blends in under clothing
- Not made with natural rubber latex
- ***Tip:*** *Measure before ordering*

Code	Size	Units	HCPCS*
4215	Up to 43⅓" length	15	A4367
4220	Up to 49" length	15	A4367

Coloplast® Urostomy Night Bag



Designed for urine output. Anti-kink tubing with universal connector; requires an adapter (included with most urostomy pouches) to connect.

- Holds up to 2 liters (2000 mL)
- 120 cm long flexible tube with anti-reflux valve
- Bag hanger (#5070) not included- order separately
- Not made with natural rubber latex

Code	Units	HCPCS*
21365	10	A4357

Coloplast® Bed Drainage Bag



Wide-lumen bag tubing specifically designed for fecal output.

- Holds up to 2 liters (2000 mL)
- Can be cut to find right length for bed drainage
- Compatible with fistula and wound management systems or high-output pouches
- Bag hanger (#5070) not included- order separately
- Not made with natural rubber latex

Code	Units	HCPCS*
14010	6	A4357

Coloplast® Bag Hanger



Reusable, sturdy, plastic drainage bag hanger secures tubing and ensures bag stability.

- Not made with natural rubber latex

Code	Size	Units	HCPCS*
5070	One Size	1	Medicare N/A Other A4421

Coloplast® Fistula Drain Port



For use with Fistula and Wound Management Systems (page 93). Drain port can be easily applied to the lid to allow for catheter insertion.

- For catheters ranging from 18FR - 22FR, the drain port is ready for use
- For catheters ranging from 24FR - 30FR, pull off the top of the drain port - DO NOT CUT
- Not made with natural rubber latex

Code	Units
14015	12

Coloplast® Transparent Fistula Lid



For use with Fistula and Wound Management Systems (page 93). Replacement lids match the transparent lid included with each fistula pouch.

- Not made with natural rubber latex

Code	Size	Units
14021	Mini (matches 14050)	10
14031	Midi (matches 14060)	10
14041	Maxi (matches 14070)	10

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Assura® Ileostomy Night Bag



Designed to collect ileostomy output generated during the post-op period. Compatible with Assura and SenSura® Click barriers.

- 2 liters (2000 mL) capacity for very high output stomas
- Clamp closure
- Comes with bag hanger
- **Filter**
- Not made with natural rubber latex

2 L — Units: 5

Code	Color	HCPCS*
Opaque		
2836		A4413
2837		A4413

Irrigation Sets

Coloplast® Irrigation Set



Durable water bag with integrated thermometer, bold print, large opening. Can be used with 1- or 2-piece ostomy systems.

- Not made with natural rubber latex

Code	Units	Includes	HCPCS*
1500	1	1- Water Bag/Regulator 1- Cone (#1110)	A4398 A4399

PDAC Product Classification List lists additional codes (A4361 and A4367) which are not included in the set.

SenSura® Mio Hospital Irrigation Set



Durable water bag with integrated thermometer, bold print, large opening.

- Not made with natural rubber latex

Code	Units	Includes	HCPCS*
1500H	1	1- Water Bag/Regulator 1- Cone (#1110) 4- SenSura® Mio 1-Piece Irrigation Sleeve, cut-to-fit with wire closure $^{3}/_{8}$ – 2 $^{1}/_{8}$" (15-55 mm) (#11750) 2- Shirt clips	A4398 A4399 A4437

PDAC Product Classification List lists additional codes (A4361, A4367, A4436) which are not included in the set.

SenSura® or Assura® Irrigation Set-Deluxe Version



Components supplied in a travel bag.

- Not made with natural rubber latex

Code	Units	Includes	HCPCS*
12830	1	1- Water Bag/Regulator 1- Cone (#1110) 1- Irrigation Faceplate (#12820) 1- Belt 2- Irrigation Sleeves (#12836) with wire closure	A4398 A4399 A4361 A4367 A4436

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Irrigation Sleeves

SenSura® Mio 1-Piece Irrigation Sleeve



1-piece irrigation sleeve with SenSura Mio adhesive is to be used by people with a colostomy who irrigate their bowels. Cut to fit stoma size.

- Elastic, skin-friendly SenSura Mio adhesive with BodyFit Technology
- Adheres to body contours creating a comfortable and secure fit
- Easy to remove without causing skin irritation
- Customizable stoma opening
- Tapered sleeve shape with wire closure
- Not made with natural rubber latex

Transparent

Code	Stoma Size	Units	HCPCS*
11750	3/8 - 2 1/8" (15-55 mm)	10	A4437

Irrigation Sleeves (without adhesive)



Irrigation sleeves with wire closure for use with SenSura® Click or Assura® barriers or faceplate. Consult chart below to order the correct transparent sleeves.

- Not made with natural rubber latex

Code	Size	Units	HCPCS*
12834		5	A4436
12835		5	A4436
12836		5	A4436

Code	Use with SenSura® Click Product Code	Units
12834	10011, 10012, 10013, 10015-10017, 11011, 11012, 11015, 11016	5
12835	10021, 10022, 10023, 10025-10027, 11021-11025	5
12836	10031, 10032, 10035, 10036, 11031, 11032, 11035, 11036	5

Irrigation Accessories

Assura® Irrigation Faceplate



For use with Brava® Belt 4215 or 4220. (Use with #12836 Irrigation Sleeve) Not compatible with Brava® Belt for SenSura® Mio.

- Not made with natural rubber latex

Code	Units	HCPCS*
12820	1	A4361

Coloplast® Stoma Cone



Soft and flexible stoma cone with tubing.

- Not made with natural rubber latex

Code	Units	HCPCS*
1110	1	A4399

SenSura® Mio

Because *every fit* matters

- All **SenSura® Mio** barriers feature **BodyFit Technology** – a unique, elastic adhesive that quickly adheres to the skin and flows into creases and folds
- **SenSura Mio** provides a secure, comfortable fit that allows your patients to feel secure and confident




Lucie
User | SenSura® Mio and Brava®

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

SenSura® Mio Click

Intuitive and **easy** to apply – with a **secure** and **leak-tight** coupling*



All it takes is a gentle press and lock to secure!

Place

1. Make sure the lock is open and **align the pouch and coupling.**

Tip: Move the lock to the 11 o'clock or 1 o'clock position so it's easy to find.



Click

2. **Apply gentle pressure** to the pouch with your hand to 'click' the coupling into place.



Tip: Start with as little pressure as possible.

Lock

3. **Squeeze the lock** until you hear a 'click'.

Tip: Audible feedback offers reassurance pouch is secure.



*Compared to prior Click coupling

NOTE: SenSura® Mio Click is **NOT** compatible with SenSura® Click or Assura® systems.

SenSura® Mio Click Flat Barrier



- Extended wear
- Flat
- Belt tabs for use with 4237, 4247
- Not made with natural rubber latex

Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
10502	⅜-1⅜" (10-35 mm)		A4409
10512	⅜-1¾" (10-45 mm)		A4410
10522	⅜-2⅛" (10-55 mm)		A4410
10531	⅜-2 9/16" (10-65 mm)		A4410
Pre-cut			
10503	¾" (20 mm)		A4409
10506	⅞" (22 mm)		A4409
10504	1" (25 mm)		A4409
10505	1⅛" (28 mm)		A4409
10513	1 3/16" (30 mm)		A4410
10514	1¼" (32 mm)		A4410
10515	1⅜" (35 mm)		A4410
10516	1½" (38 mm)		A4410
10523	1 9/16" (40 mm)		A4410
10524	1¾" (45 mm)		A4410
10525	2" (50 mm)		A4410

SenSura® Mio Click Convex Light Barrier



- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Not made with natural rubber latex

Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
16901	⅜ - ¾" (10-20 mm)		A4407
16911	⅝ - 1 3/16" (15-30 mm)		A4408
16921	⅝ - 1 9/16" (15-40 mm)		A4408
16931	⅝ - 2 1/16" (15-53 mm)		A4408
Pre-cut			
16903	11/16" (18 mm)		A4407
16913	13/16" (21 mm)		A4408
16914	1" (25 mm)		A4408
16915	1⅛" (28 mm)		A4408
16923	1¼" (31 mm)		A4408
16924	1⅜" (35 mm)		A4408
16933	1½" (38 mm)		A4408
16934	1⅝" (41 mm)		A4408

Convex Light

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

SenSura® Mio Click Convex Deep Barrier

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Not made with natural rubber latex



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
16941	3/8 - 3/4" (10-20 mm)		A4407
16951	5/8 - 1 3/16" (15-30 mm)		A4408
16961	5/8 - 1 3/16" (15-40 mm)		A4408
16971	5/8 - 2 1/16" (15-53 mm)		A4408
Pre-cut			
16942	11/16" (18 mm)		A4407
16953	13/16" (21 mm)		A4408
16954	1" (25 mm)		A4408
16955	1 1/8" (28 mm)		A4408
16963	1 1/4" (31 mm)		A4408
16964	1 3/8" (35 mm)		A4408
16973	1 1/2" (38 mm)		A4408
16974	1 5/8" (41 mm)		A4408

Convex Deep

NOTE: SenSura® Mio Click is **NOT** compatible with SenSura® Click or Assura® systems.

SenSura® Mio Click Convex Flip Barrier

- Extended wear
- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- Belt tabs for use with 4237, 4247
- Not made with natural rubber latex



**SenSura Mio Convex Flip is available in multiple barrier sizes to conform to different sizes of outward areas. Start with the smallest barrier size that fits the stoma. For larger outward areas, size up as needed.
- S: Max cut of 1 9/16" (40 mm)
- M: Max cut of 2" (50 mm)
- L: Max cut of 2 3/8" (60 mm)
- XL: Max cut of 2 5/8" (66 mm)

Units: 5

Code	Stoma Size	Barrier size**	Coupling size	HCPCS*
Cut-to-fit				
18500	3/8 - 1 3/8" (10-35 mm)	S		A4407
18510	3/8 - 1 3/4" (10-45 mm)	M		A4408
18520	3/8 - 2 1/8" (10-55 mm)	L		A4408
18530	3/8 - 2 9/16" (10-66 mm)	XL		A4408

Code	Stoma Size	Barrier size**	Coupling size	HCPCS*
Pre-cut				
18501	3/4" (20 mm)	S		A4407
18511	1" (25 mm)	M		A4407
18512	1 3/16" (30 mm)	M		A4407
18513	1 3/8" (35 mm)	M		A4407
18521	1 3/16" (30 mm)	L		A4407
18522	1 3/8" (35 mm)	L		A4407
18523	1 9/16" (40 mm)	L		A4407
18524	1 3/4" (45 mm)	L		A4407
18531	1 3/4" (45 mm)	XL		A4408
18532	2" (50 mm)	XL		A4408

Indicates Convex Flip

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

SenSura® Mio Click MAXI Drainable Pouch



Improved coupling is as secure, but easier to apply - with less pressure^

- EasiClose™ WIDE outlet
- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile
- Opaque includes inspection window
- Not made with natural rubber latex

11" (28 cm) / 655 mL — Units: 10

Code	Code	Code	Color	HCPCS*
Transparent	Opaque gray	Opaque black		
11463	**11462**	**19759**		A4427
11473	**11472**	**19761**		A4427
11483	**11482**	**19762**		A4427
11487	**11486**	**19763**		A4427

SenSura® Mio Click MAXI Drainable Pouch



Improved coupling is as secure, but easier to apply - with less pressure^

- EasiClose™ WIDE outlet
- **Without filter**
- Neutral gray, water-resistant textile
- Not made with natural rubber latex

11" (28 cm) / 655 mL — Units: 10

Code	Color	HCPCS*
Transparent		
11474		A4426
11484		A4426
11488		A4426

SenSura® Mio Click MIDI Drainable Pouch



Improved coupling is as secure, but easier to apply - with less pressure^


- EasiClose™ WIDE outlet
- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile
- Opaque includes inspection window
- Not made with natural rubber latex

10¼" (26 cm) / 485 mL — Units: 10

Code	Code	Color	HCPCS*
Opaque gray	Opaque black		
11461	**19758**		A4427
11471	**19760**		A4427

^Compared to prior Click coupling.
Coloplast, Data-on-file (09/2023), NCT05652010, n=39.
Simulated use by users and HCPs, Coloplast Data-on-file (03/2024), n=30; 97% [95%CI: 83%-100%]

SenSura® Mio Click High Output Pouch features



Transparent
- Printed scale for visual indication of output volume

Opaque
- Inspection window
- Option to tuck up outlet for discretion

SenSura® Mio Click High Output Pouch with TAP Outlet



Improved coupling is as secure, but easier to apply - with less pressure^

- **Full-circle filter designed for ileostomy output**
- Neutral gray, water-resistant textile
- Not considered for use with fully formed stools. Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast® Bed Drainage Bag #14010
- Not made with natural rubber latex

12½" (32 cm) / 1080 mL — Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque gray		
18610	**18600**		A4413
18611	**18601**		A4413
18612	**18602**		A4413
18613	**18603**		A4413

SenSura® Mio Click High Output Pouch with SOFT Outlet



Improved coupling is as secure, but easier to apply - with less pressure^

- **Full-circle filter designed for ileostomy output**
- Neutral gray, water-resistant textile
- For thin, liquid, high-output stool
- Soft outlet provides effective drainage; connects to a universal bedside drainage bag or Coloplast® Bed Drainage Bag #14010
- Not made with natural rubber latex

12½" (32 cm) / 1080 mL — Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque gray		
18640	**18630**		A4413
18641	**18631**		A4413
18642	**18632**		A4413

^Compared to prior Click coupling.
Coloplast, Data-on-file (09/2023), NCT05652010, n=39.
Simulated use by users and HCPs, Coloplast Data-on-file (03/2024), n=30; 97% [95%CI: 83%-100%]

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

SenSura® Mio Click Drainable Pouch with TAP Outlet



Improved coupling is as secure, but easier to apply - with less pressure^

- **Full-circle filter designed for ileostomy output**
- Neutral gray, water-resistant textile
- Opaque includes inspection window
- Not considered for use with fully formed stools. Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast® Bed Drainage Bag #14010
- Not made with natural rubber latex

9" (23 cm) / 650 mL — Units: 10

Code	Color	HCPCS*
Opaque gray		
18620		A4427
18621		A4427
18622		A4427

SenSura® Mio Click MAXI Closed Pouch



Improved coupling is as secure, but easier to apply - with less pressure^

- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile
- Opaque includes inspection window
- Not made with natural rubber latex

8¼" (20½ cm) / 570 mL — Units: 30

Code	Code	Code	Color	HCPCS*
Transparent	Opaque gray	Opaque black		
N/A	**11402**	**19727**		A4423
N/A	**11413**	**19730**		A4423
N/A	**11422**	**19732**		A4423
11426	**11425**	**19733**		A4423

SenSura® Mio Click MIDI Closed Pouch



Improved coupling is as secure, but easier to apply - with less pressure^

- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile
- Opaque includes inspection window
- Not made with natural rubber latex

7¼" (18½ cm) / 275 mL — Units: 30


Code	Code	Color	HCPCS*
Opaque gray	Opaque black		
11401	**19726**		A4423
11412	**19729**		A4423
11421	**19731**		A4423

^Compared to prior Click coupling.
Coloplast, Data-on-file (09/2023), NCT05652010, n=39.
Simulated use by users and HCPs, Coloplast Data-on-file (03/2024), n=30; 97% [95%CI: 83%-100%]

2-Piece SenSura® Mio Click

Urostomy

SenSura® Mio Urostomy Pouch with Twist Outlet



Wide anti-reflux valves support varying ways of managing stents:

- Threading the stents into the pouch via the anti-reflux valves
- Rolling the stents into a loop at the top of the pouch
- Trimming the stents

SenSura® Mio Click MAXI Urostomy Pouch with Twist Outlet



Improved coupling is as secure, but easier to apply - with less pressure^

- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Twist outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Neutral gray, water-resistant textile
- Not made with natural rubber latex

10½" (26½ cm) / 480 mL — Units: 10 pouches, 5 adapters

Code	Code	Color	HCPCS*
Transparent	Opaque gray		
114920	**114910**		A4434
114980	**114970**		A4434
115000	**114990**		A4434

Urostomy Adapter



- Urostomy adapter - connect the new SenSura® Mio urostomy pouch's outlet to any universal bedside drainage bag
- Easy clip-on mechanism - secures to the outlet with a gentle pinch enabling connection to bedside drainage bag

Code	Units	HCPCS*
18790*	5 adapters	NA

*Purchase available for hospital / acute settings only.
Patients who need additional adapters can look to request for samples by calling Coloplast Care 1-877-858-2656

^Compared to prior Click coupling

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Single Use Kits

SenSura® Mio Click Flat Barrier + Drainable Pouch Single Use Kit

Improved coupling is as secure, but easier to apply - with less pressure^

- Extended wear + pouch
- Instructions for use
- Stoma measuring guide



Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
13651	⅜ - 1¾" (10-45 mm)	
Barrier #10512 Flat, Pouch #11473 MAXI / EasiClose™ WIDE outlet / filter		
21008	⅜ - 2⅛" (10-55 mm)	
Barrier #10522 Flat, Pouch #11483 MAXI / EasiClose™ WIDE outlet / filter		
13652	⅜ - 2½" (10-65 mm)	
Barrier #10531 Flat, Pouch #11487 MAXI / EasiClose™ WIDE outlet / filter		

SenSura® Mio Click Convex Light Barrier + Drainable Pouch Single Use Kit

Improved coupling is as secure, but easier to apply - with less pressure^

- Extended wear + pouch
- Instructions for use
- Stoma measuring guide



Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
13653	⅝ - 1⅛" (15-30 mm)	
Barrier #16911 Convex Light, Pouch #11473 MAXI / EasiClose™ WIDE outlet / filter		
13654	⅝ - 1⅝" (15-40 mm)	
Barrier #16921 Convex Light, Pouch #11483 MAXI / EasiClose™ WIDE outlet / filter		
13655	⅝ - 2 1/16" (15-53 mm)	
Barrier #16931 Convex Light, Pouch #11487 MAXI / EasiClose™ WIDE outlet / filter		

Convex Light

SenSura® Mio Click Flat Barrier + MAXI Urostomy Pouch with Twist Outlet Single Use Kit

Improved coupling is as secure, but easier to apply - with less pressure^

- Extended wear + pouch
- Instructions for use
- Stoma measuring guide
- Urostomy adapter



Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
13656	⅜ - 1¾" (10-45 mm)	
Barrier #10512 Flat, Pouch #114980 Twist outlet / multi-chamber		
13657	⅜ - 2¼" (10-55 mm)	
Barrier #10522 Flat, Pouch #115000 Twist outlet / multi-chamber		

SenSura® Mio Click Convex Light Barrier + MAXI Urostomy Pouch with Twist Outlet Single Use Kit

Improved coupling is as secure, but easier to apply - with less pressure^

- Extended wear + pouch
- Instructions for use
- Stoma measuring guide
- Urostomy adapter



Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
13658	⅝ - 1 3/16" (15-30 mm)	
Barrier #16911 Convex Light, Pouch #114980 Twist outlet / multi-chamber		
13659	⅝ - 1 9/16" (15-40 mm)	
Barrier #16921 Convex Light, Pouch #115000 Twist outlet / multi-chamber		

Convex Light

^Compared to prior Click coupling

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

NOTE: SenSura® Mio Click is NOT compatible with SenSura® Click or Assura® systems. SenSura® Click and Assura® are compatible.

SenSura® Click Flat Barrier



- Standard wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex

Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
10011	⅜ - 1⅜" (10-35 mm)		A4414
10021	⅜ - 1¾" (10-45 mm)		A4415
10031	⅜ - 2⅛" (10-55 mm)		A4415
10041	⅜ - 2½" (10-65 mm)		A4415
Pre-cut			
10012	13⁄16" (20 mm)		A4414
10013	1" (25 mm)		A4414
10022	1 3⁄16" (30 mm)		A4415
10023	1⅜" (35 mm)		A4415
10032	1 9⁄16" (40 mm)		A4415

SenSura® Click Xpro Flat Barrier



- Extended wear
- Belt tabs
- Not made with natural rubber latex

Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
10015	⅜ - 1⅜" (10-35 mm)		A4409
10025	⅜ - 1¾" (10-45 mm)		A4410
10035	⅜ - 2¼" (10-55 mm)		A4410
10045	⅜ - 2½" (10-65 mm)		A4410
Pre-cut			
10016	13⁄16" (20 mm)		A4409
10017	1" (25 mm)		A4409
10026	1 3⁄16" (30 mm)		A4410
10027	1⅜" (35 mm)		A4410
10036	1 9⁄16" (40 mm)		A4410

SenSura® Click Convex Light Barrier



- Standard wear
- Belt tabs
- Not made with natural rubber latex

Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
11011	⅝ - ⅞" (15-22 mm)		A4373
11021	⅝ - 1 5⁄16" (15-33 mm)		A4373
11031	⅝ - 1 11⁄16" (15-43 mm)		A4373
11041	⅝ - 2 1⁄16" (15-53 mm)		A4373
Pre-cut			
11012	⅞" (22 mm)		A4373
11022	1" (25 mm)		A4373
11023	1⅛" (28 mm)		A4373
11024	1¼" (31 mm)		A4373
11032	1⅜" (35 mm)		A4373

SenSura® Click Xpro Convex Light Barrier



- Extended wear
- Belt tabs
- Not made with natural rubber latex

Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
11015	⅝ - ⅞" (15-22 mm)		A4407
11025	⅝ - 1 5⁄16" (15-33 mm)		A4408
11035	⅝ - 1 11⁄16" (15-43 mm)		A4408
11045	⅝ - 2 1⁄16" (15-53 mm)		A4408
Pre-cut			
11016	⅞" (22 mm)		A4407
11026	1" (25 mm)		A4408
11027	1⅛" (28 mm)		A4408
11028	1¼" (31 mm)		A4408
11036	1⅜" (35 mm)		A4408

Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

SenSura® Click MAXI Drainable Pouch



- EasiClose™ WIDE outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11½" (29½ cm) / 655 mL Units: 20

Code	Code	Color	HCPCS*
Transparent	Opaque		
11194	**11124**		A4427
11195	**11125**		A4427
11196	**11126**		A4427
11197	**11127**		A4427**

SenSura® Click MIDI Drainable Pouch



- EasiClose™ WIDE outlet
- **Filter**
- Soft cloth on front and back
- Not made with natural rubber latex

10½" (26½ cm) / 485 mL Units: 20

Code	Color	HCPCS*
Opaque		
11114		A4427
11115		A4427

NOTE: SenSura® Click is compatible with Assura® but NOT with SenSura® Mio Click.

SenSura® Click MAXI Drainable Pouch



- EasiClose™ WIDE outlet
- **Without filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11½" (29½ cm) / 655 mL Units: 20

Code	Code	Color	HCPCS*
Transparent	Opaque		
11185	**11135**		A4426
11186	**11136**		A4426

SenSura® Click Magnum High Output Drainable Pouch



- Not considered for use with fully formed stools. Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast® Bed Drainage Bag #14010
- **Filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

12½" (32 cm) / 1160 mL Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
19050	**19040**		A4413**
19051	**19041**		A4413**
19052	**19042**		A4413**

For Assura® Ileostomy Night Bag, see page 24.

SenSura® Click MAXI Closed Pouch



- **Filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

8½" (21 cm) / 570 mL Units: 30

Code	Code	Color	HCPCS*
Transparent	Opaque		
NA	**10164**		A4423
NA	**10165**		A4423
NA	**10166**		A4423
10187	**10167**		A4423**

SenSura® Click MIDI Closed Pouch



- **Filter**
- Soft cloth on front and back
- Not made with natural rubber latex

7" (18 cm) / 275 mL Units: 30

Code	Color	HCPCS*
Opaque		
10154		A4423
10155		A4423
10156		A4423

For Assura® Stoma Caps, see page 50.

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

**Product has not been through PDAC verification

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Urostomy

SenSura® Click MAXI Urostomy Pouch



- Multi-chambers minimize sloshing
- Wider anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

10⅜" (26 cm) / 480 mL — Units: 10 pouches, 5 adapters

Code	Code	Color	HCPCS*
Transparent	Opaque		
11854	**11844**		A4433
11855	**11845**		A4433
11856	**11846**		A4433

SenSura® Click MIDI Urostomy Pouch



- Multi-chambers minimize sloshing
- Wider anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

9½" (24 cm) / 430 mL — Units: 10 pouches, 5 adapters

Code	Code	Color	HCPCS*
Transparent	Opaque		
11851	**11841**		A4433

Single Use Kits

SenSura® Click Flat Barrier + Drainable Pouch Single Use Kit



- Standard wear + pouch
- Instructions for use
- Stoma measuring guide

Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
11009	⅜ - 2¼" (10-55 mm)	

Barrier #10031 Flat, Pouch #11196 MAXI / EasiClose™ WIDE outlet / filter

SenSura® Click Convex Light Barrier + Drainable Pouch Single Use Kit



- Standard wear + pouch
- Instructions for use
- Stoma measuring guide

Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
11010	⅝ - 1 11/16" (15-43 mm)	

Barrier #11031 Convex Light, Pouch #11196 MAXI / EasiClose™ WIDE outlet / filter

SenSura® Click Xpro Flat, Extended Wear Barrier + Urostomy Pouch Single Use Kit



- Urostomy adapter
- Instructions for use
- Stoma measuring guide

Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
11004	⅜ - 2¼" (10-55 mm)	

Barrier #10035 Flat, Pouch #11856 MAXI / soft outlet / multi-chamber

Urostomy Night Bag and Bag Hanger

See Supporting Products on page 23 for more information.



Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

NOTE: Assura® is compatible with SenSura® Click , but **NOT** with SenSura® Mio Click.

Assura® Flat Barrier



- Standard wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex

Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
2881	⅜ - 1⅜" (10-35 mm)		A4414
2882	⅜ - 1¾" (10-45 mm)		A4414
2883	⅜ - 2¼" (10-55 mm)		A4414
Pre-cut			
12841	13/16" (20 mm)		A4414
12842	1" (25 mm)		A4414
12844	1 3/16" (30 mm)		A4414
12845	1⅜" (35 mm)		A4414
12847	1 3/16" (40 mm)		A4414
12848	1¾" (45 mm)		A4414

Assura® Flat Barrier



- Extended wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex

Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
2831	⅜ - 1⅜" (10-35 mm)		A4409
2832	⅜ - 1¾" (10-45 mm)		A4409
2833	⅜ - 2¼" (10-55 mm)		A4409
Pre-cut			
14250	⅝" (15 mm)		A4409
14251	13/16" (20 mm)		A4409
14252	1" (25 mm)		A4409
14254	1 3/16" (30 mm)		A4409
14255	1⅜" (35 mm)		A4409
14257	1 3/16" (40 mm)		A4409
14258	1¾" (45 mm)		A4409

Assura® Convex Light Barrier



- Standard wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex

Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
14261	5/8 - 15/16" (15-23 mm)		A4373
14262	5/8 - 1 5/16" (15-33 mm)		A4373
14263	5/8 - 1 11/16" (15-43 mm)		A4373
Pre-cut			
14271	¾" (18 mm)		A4373
14272	⅞" (22 mm)		A4373
14273	1" (25 mm)		A4373
14274	1⅛" (28 mm)		A4373
14275	1¼" (31 mm)		A4373
14276	1⅜" (35 mm)		A4373
14277	1½" (38 mm)		A4373
14278	1⅝" (41 mm)		A4373

Assura® Convex Light Barrier



- Extended wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex

Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
14281	5/8 - 15/16" (15-23 mm)		A4407
14282	5/8 - 1 5/16" (15-33 mm)		A4407
14283	5/8 - 1 11/16" (15-43 mm)		A4407
Pre-cut			
14291	¾" (18 mm)		A4407
14292	⅞" (22 mm)		A4407
14293	1" (25 mm)		A4407
14294	1⅛" (28 mm)		A4407
14295	1¼" (31 mm)		A4407
14296	1⅜" (35 mm)		A4407
14297	1½" (38 mm)		A4407
14298	1⅝" (41 mm)		A4407

Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

NOTE: Assura® is compatible with SenSura® Click, but **NOT** with SenSura® Mio Click.

Assura® Deep Convex Barrier



- Standard wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex

Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
12716	⅝ - 1 5/16" (15-33 mm)		A4373
12719	⅝ - 1 11/16" (15-43 mm)		A4373
Pre-cut			
12701	⅝" (15 mm)		A4373
12702	¾" (18 mm)		A4373
12703	⅞" (22 mm)		A4373
12704	1" (25 mm)		A4373
12705	1⅛" (28 mm)		A4373
12706	1¼" (31 mm)		A4373
12707	1⅜" (35 mm)		A4373
12708	1½" (38 mm)		A4373
12709	1⅝" (41 mm)		A4373

Assura® Deep Convex Barrier



- Extended wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex

Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
14243	5/8 - 15/16" (15-23 mm)		A4407
14246	5/8 - 1 5/16" (15-33 mm)		A4407
14249	5/8 - 1 11/16" (15-43 mm)		A4407
Pre-cut			
14232	¾" (18 mm)		A4407
14233	7/8" (22 mm)		A4407
14234	1" (25 mm)		A4407
14235	1 1/8" (28 mm)		A4407
14236	1¼" (31 mm)		A4407
14237	1⅜" (35 mm)		A4407
14238	1½" (38 mm)		A4407
14239	1 5/8" (41 mm)		A4407

Indicates convexity

Assura® MAXI Drainable Pouch



- EasiClose™ WIDE outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11½" (29½ cm) / 655 mL Units: 10

Code		Color	HCPCS*
Transparent	Opaque		
15974	**15984**		A4427
15975	**15985**		A4427
15976	**15986**		A4427

Assura® MAXI Drainable Pouch



- EasiClose™ WIDE outlet
- **Without filter**
- Soft cloth on front and back
- Not made with natural rubber latex

11½" (29½ cm) / 655 mL Units: 10

Code	Color	HCPCS*
Opaque		
15982		A4426
15983		A4426

Assura® MIDI Drainable Pouch



- EasiClose™ WIDE outlet
- **Filter**
- Soft cloth on front and back
- Not made with natural rubber latex

10½" (26½ cm) / 485 mL Units: 10

Code	Color	HCPCS*
Opaque		
15954		A4427
15955		A4427
15956		A4427

Assura® MAXI Drainable Pouch



- EasiClose™ outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11½" (29½ cm) / 655 mL Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
13974	**13984**		A4427
13975	**13985**		A4427
13976	**13986**		A4427

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

**Product has not been through PDAC verification

NOTE: Assura® is compatible with SenSura® Click, but **NOT** with SenSura® Mio Click.

Assura® MIDI Drainable Pouch



- EasiClose™ outlet
- **Filter**
- Soft cloth on front and back
- Graphic pattern
- Not made with natural rubber latex

10½" (26½ cm) / 485 mL — Units: 10

Code	Color	HCPCS*
Opaque		
13964		A4427
13965		A4427
13966		A4427

Assura® MINI Drainable Pouch



- EasiClose™ outlet
- **Filter**
- Soft cloth on front and back
- Not made with natural rubber latex

9¼" (23½ cm) / 280 mL — Units: 10

Code	Color	HCPCS*
Opaque		
13924		A4427
13925		A4427

Assura® MAXI Drainable Pouch



- Clamp closure
- **Without filter**
- Soft cloth on front and back
- Not made with natural rubber latex

11¼" (28½ cm) / 760 mL — Units: 10

Code	Color	HCPCS*
Opaque		
12577		A4426
12578		A4426
12579		A4426

Assura® MIDI Drainable Pouch



- Clamp closure
- **Without filter**
- Soft cloth on front and back
- Not made with natural rubber latex

9¾" (24½ cm) / 400 mL — Units: 10

Code	Color	HCPCS*
Opaque		
12527		A4426
12528		A4426
12529		A4426

Assura® High Output Drainable Pouch



- Soft outlet can connect to Coloplast® Bed Drainage Bag #14010, or a drain bag with a standard connection
- **Filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

12½" (32 cm) / 750 mL — Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
8114	**NA**		A4413
8115	**2846**		A4413
8116	**2847**		A4413

For urostomy night bag and bag hanger, see page 23.

Assura® Ileostomy Night Bag



Designed to collect ileostomy output generated during the post-op period. Compatible with Assura and SenSura® Click barriers.

- 2 liters (2000 mL) capacity for very high output stomas
- Clamp closure
- Comes with bag hanger
- **Filter**
- Not made with natural rubber latex

2 L — Units: 5

Code	Color	HCPCS*
Opaque		
2836		A4413
2837		A4413

Assura® MAXI Closed Pouch



- **Filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

8½" (21 cm) / 570 mL — Units: 30

Code	Code	Color	HCPCS*
Transparent	Opaque		
12374	**12384**		A4423
12375	**12385**		A4423
12376	**12386**		A4423

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Assura® MIDI Closed Pouch



- **Filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

7" (18 cm) / 275 mL — Units: 30

Code	Code	Color	HCPCS*
Transparent	Opaque		
12344	**12354**		A4423
12345	**12355**		A4423
12346	**12356**		A4423

Ileostomy / Colostomy

Assura® Stoma Cap



- **Filter**
- Small size is optimal for intimate moments or colostomates who irrigate
- Not made with natural rubber latex

Units: 30

Code	Color	HCPCS*
Opaque		
2801		A5055
2802		A5055
2803		A5055**

Urostomy

Assura® Urostomy Micro-Pouch



- Miniature size prevents twisting
- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for patients in a wheelchair
- Soft cloth on front
- Not made with natural rubber latex

6" (15 cm) — Units: 10

Code	Color	HCPCS*
Transparent		
14204		A4433
14205		A4433
14206		A4433

For nighttime, use a urostomy night bag found on page 23.

**Product has not been through PDAC verification

Assura® Urostomy Minicap



- Ultra-absorbent material inside
- Small size is optimal for intimate use
- Not made with natural rubber latex

100 mL — Units: 30

Code	Color	HCPCS*
Opaque		
2807		A5055
2808		A5055

Assura® MAXI Urostomy Pouch



- Multi-chambers minimize sloshing
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

10½" (26½ cm) / 488 mL — Units: 10 pouches, 5 adapters

Code	Code	Color	HCPCS*
Transparent	Opaque		
14227	**14224**		A4433
14228	**14225**		A4433
14229	**14226**		A4433

Assura® MAXI Urostomy Pouch



- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

10¾" (27¼ cm) / 375 mL — Units: 10 pouches, 5 adapters

Code	Code	Color	HCPCS*
Transparent	Opaque		
1757	**1754**		A4433
1758	**1755**		A4433
1759	**1756**		A4433

Assura® MIDI Urostomy Pouch



- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Soft cloth on back
- Not made with natural rubber latex

9½" (24 cm) / 150 mL — Units: 10 pouches, 5 adapters

Code	Color	HCPCS*
Transparent		
1751		A4433
1752		A4433

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Single Use Kits

Assura® 2-piece Flat Barrier+ MAXI Drainable Pouch Single Use Kit
- Standard wear + pouch
- Instructions for use
- Stoma measuring guide



Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
11002	⅜ - 2¼" (10-55 mm)	

Barrier #2883 Flat, Pouch #15976 MAXI / EasiClose™ WIDE outlet / filter

Assura® 2-piece Convex Light Barrier + MAXI Drainable Pouch Single Use Kit
- Standard wear + pouch
- Instructions for use
- Stoma measuring guide



Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
11001	⅝ - 1 11/16" (15-43 mm)	

Barrier #14263 Convex Light, Pouch #15976 MAXI / EasiClose™ WIDE outlet / filter

Indicates convexity

SenSura®Mio

Patients deserve proven leakage reduction[1] from the start

A multi-site US acute care study showed

SenSura® Mio users had 41% lower risk of leakage

compared to users of Hollister® CeraPlus™ pouching systems[1*]



1. Based on a retrospective study conducted at 2 affiliated university hospitals; n=214 patient charts with surgery dates ranging from February 2022 to June 2023. Patient events were charted by the WOC nurse at each pouch change; review captured up to 5 pouch changes or 2 weeks of data. Statistical analysis controlled for key differences between sites and found predicted probability of leakage for SenSura Mio was 15.1% and Hollister® CeraPlus™ pouching system was 25.6%, corresponding to a 40.8% lower risk (p=0.011). Zamarripa, C. et al. Pouching System Leakage and Peristomal Skin Complications Following Ostomy Surgery in the Immediate Postoperative Period: A Retrospective Review. *J Wound Ostomy Continence Nurs.* 2024;51(6):478-483.

™*All trademarks are the property of their respective owners

SenSura® Mio Flex

All the benefits of BodyFit Technology with a 2-piece adhesive coupling system

- The Flex coupling is a comfortable, flexible, low-profile 2-piece system
- The floating flange allows the pouch to be applied to the barrier with minimal pressure on the abdomen





NOTE: SenSura® Mio Flex and SenSura® Flex **are compatible** with one another

SenSura® Mio Flex Flat Barrier

- Extended wear
- Not made with natural rubber latex



Units: 5
18390† - 10 units

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
10551	⅜ - 1 5/16" (10-33 mm)		A4410
10561	⅜ - 1⅞" (10-48 mm)		A4410
10571	⅜ - 2 11/16" (10-68 mm)		A4410
18390†	⅜ - 3½" (10-88 mm)		A4410
Pre-cut			
10562	1" (25 mm)		A4410
10566	1⅛" (28 mm)		A4410
10563	1 3/16" (30 mm)		A4410
10567	1¼" (32 mm)		A4410
10564	1⅜" (35 mm)		A4410
10568	1½" (38 mm)		A4410
10565	1 9/16" (40 mm)		A4410
10572	1¾" (45 mm)		A4410
10573	2" (50 mm)		A4410

SenSura® Mio Flex Convex Soft Barrier

Gentle support: highly compressible, highly flexible

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Not made with natural rubber latex



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
16880	⅜ - ¾" (10-20 mm)		A4408
16882	⅝ - 1 1/8" (15-30 mm)		A4408
16886	⅝ - 1⅝" (15-40 mm)		A4408
16889	⅝ - 2 1/16" (15-53 mm)		A4408
Pre-cut			
16881	11/16" (18 mm)		A4408
16883	13/16" (21 mm)		A4408
16884	1" (25 mm)		A4408
16885	1⅛" (28 mm)		A4408
16887	1¼" (31 mm)		A4408
16888	1⅜" (35 mm)		A4408
16890	1½" (38 mm)		A4408
16891	1⅝" (41 mm)		A4408

Convex Soft

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

NOTE: SenSura® Mio Flex and SenSura® Flex **are compatible** with one another

SenSura® Mio Flex Convex Light Barrier



Light support: moderately compressible, highly flexible
- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Not made with natural rubber latex

Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
16471	⅜ - ¾" (10-20 mm)		A4408
16481	⅝ - 1 3/16" (15-30 mm)		A4408
16483	⅝ - 1 9/16" (15-40 mm)		A4408
16491	⅝ - 2 1/16" (15-53 mm)		A4408
Pre-cut			
16472	11/16" (18 mm)		A4408
16484	13/16" (21 mm)		A4408
16485	1" (25 mm)		A4408
16486	1⅛" (28 mm)		A4408
16487	1¼" (31 mm)		A4408
16488	1⅜" (35 mm)		A4408
16493	1½" (38 mm)		A4408
16494	1⅝" (41 mm)		A4408

SenSura® Mio Flex Convex Deep Barrier



Deepest on the market! Less compressible yet still highly flexible
- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Not made with natural rubber latex

Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
18540	⅜ - ¾" (10-20 mm)		A4408**
18542	⅝ - 1 3/16" (15-30 mm)		A4408**
18546	⅝ - 1⅝" (15-40 mm)		A4408**
18549	⅝ - 2 1/16" (15-53 mm)		A4408**
Pre-cut			
18541	11/16" (18 mm)		A4408**
18543	13/16" (21 mm)		A4408**
18544	1" (25 mm)		A4408**
18545	1⅛" (28 mm)		A4408**
18547	1¼" (31 mm)		A4408**
18548	1⅜" (35 mm)		A4408**
18550	1½" (38 mm)		A4408**
18551	1⅝" (41 mm)		A4408**

Convex Light Convex Deep

SenSura® Mio Flex MAXI Drainable Pouch



- EasiClose™ WIDE outlet
- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile with inspection window
- Not made with natural rubber latex

11" (28 cm) / 655 mL Units: 10

Code	Code	Code	Color	HCPCS*
Transparent	Opaque gray	Opaque black		
NA	**12262**	**19792**		A4425
12273	**12272**	**19794**		A4425
12283	**12282**	**19795**		A4425
18389	**NA**	**NA**		A4425

SenSura® Mio Flex MAXI Drainable Pouch



- EasiClose™ WIDE outlet
- Without filter
- Neutral gray, water-resistant textile
- Not made with natural rubber latex

11" (28 cm) / 655 mL Units: 10

Code	Color	HCPCS*
Transparent		
12274		A5063
12284		A5063

SenSura® Mio Flex MIDI Drainable Pouch



- EasiClose™ WIDE outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile with inspection window
- Not made with natural rubber latex

10¼" (26 cm) / 485 mL Units: 10

Code	Color	HCPCS*
Opaque gray		
12281		A4425


*For important information regarding the use of HCPCS codes listed and code definitions see page 124
**Product has not been through PDAC verfication

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

SenSura® Mio Flex High Output Pouch features

Transparent
- Printed scale for visual indication of output volume

Opaque
- Inspection window
- Option to tuck up outlet for discretion



SenSura® Mio Flex High Output Pouch with TAP Outlet



- **Full-circle filter designed for ileostomy output**
- Not considered for use with fully formed stools. Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast® Bed Drainage Bag #14010
- Not made with natural rubber latex

12½" (32 cm) / 1080 mL Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque gray		
18655	**18650**		A4413
18656	**18651**		A4413

NOTE: SenSura® Mio Flex and SenSura® Flex **are compatible** with one another

SenSura® Mio Flex High Output Pouch with SOFT Outlet



- **Full-circle filter designed for ileostomy output**
- For thin, liquid, high-output stool
- Soft outlet provides effective drainage; connects to a connects to universal fecal collection bag and Coloplast® Bed Drainage Bag #14010
- Not made with natural rubber latex

12½" (32 cm) / 1080 mL Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque gray		
18657	**18652**		A4413
18658	**18653**		A4413

SenSura® Mio Flex MAXI Closed Pouch



- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile with inspection window
- Not made with natural rubber latex

8¼" (20½ cm) / 570 mL Units: 30


Code	Code	Code	Color	HCPCS*
Transparent	Opaque gray	Opaque black		
12203	**12202**	**19769**		A4419
12216	**12213**	**19772**		A4419
12223	**12222**	**19773**		A4419

Urostomy

SenSura® Mio Urostomy Pouch with Twist Outlet

Wide anti-reflux valves support varying ways of managing stents:

- Threading the stents into the pouch via the anti-reflux valves
- Rolling the stents into a loop at the top of the pouch
- Trimming the stents



SenSura® Mio Flex MAXI Urostomy Pouch with Twist Outlet



- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Twist outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Neutral gray, water-resistant textile
- Not made with natural rubber latex

10½" (27 cm) / 400 mL Units: 10 pouches, 5 adapters

Code	Code	Color	HCPCS*
Transparent	Opaque gray		
122920	**122910**		A4432
122980	**122970**		A4432

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

NOTE: SenSura® Mio Flex and SenSura® Flex **are compatible** with one another

SenSura® Mio Flex MIDI Urostomy Pouch with Twist Outlet



- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Twist outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Neutral gray, water-resistant textile
- Not made with natural rubber latex

9½" (24 cm) / 330 mL — Units: 10 pouches, 5 adapters

Code	Color	HCPCS*
Opaque gray		
122960		A5073

Single Use Kits

SenSura® Mio Flex Flat Barrier + MAXI Drainable Pouch Single Use Kit



- Extended wear barrier + pouch
- Instructions for use
- Stoma measuring guide

Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
21009	⅜ - 1⅞" (10-48 mm)	
Barrier #10561 Flat, Pouch #12273 EasiClose™ WIDE outlet / filter		
21010	⅜ - 2¾" (10-68 mm)	
Barrier #10571 Flat, Pouch #12283 EasiClose™ WIDE outlet / filter		
13669	⅜ - 3½" (10-88 mm)	
Barrier #18390 Flat, Pouch #18389 EasiClose™ WIDE outlet / filter		
21012	⅜ - 1⅞" (10-48 mm)	
Barrier #10561 Flat, Pouch #12274 EasiClose™ WIDE outlet / **without filter**		
21013	⅜ - 2¾" (10-68 mm)	
Barrier #10571 Flat, Pouch #12284 EasiClose™ WIDE outlet / **without filter**		

SenSura® Mio Flex Convex Soft Barrier + MAXI Drainable Pouch Single Use Kit



- Extended wear barrier + pouch
- Instructions for use
- Stoma measuring guide

Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
13666	⅝ - 1⅝" (15-40 mm)	
Barrier #16886 Convex Soft, Pouch #12273 EasiClose™ WIDE outlet / filter		
13667	⅝ - 2 1/16" (15-53 mm)	
Barrier #16889 Convex Soft, Pouch #12283 EasiClose™ WIDE outlet / filter		

SenSura® Mio Flex Convex Light Barrier + MAXI Drainable Pouch Single Use Kit



- Extended wear barrier + pouch
- Instructions for use
- Stoma measuring guide

Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
13660	⅝ - 1 9/16" (15-40 mm)	
Barrier #16483 Convex Light, Pouch #12273 EasiClose™ WIDE outlet / filter		
13661	⅝ - 2 1/16" (15-53 mm)	
Barrier #16491 Convex Light, Pouch #12283 EasiClose™ WIDE outlet / filter		

SenSura® Mio Flex Flat Barrier + Urostomy Pouch with Twist Outlet Single Use Kit



- Extended wear barrier + pouch
- Urostomy adapter
- Instructions for use
- Stoma measuring guide

Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
110072	⅜ - 1⅞" (10-48 mm)	
Barrier #10561 Flat, Pouch #122980 Twist outlet / multi-chamber		

SenSura® Mio Flex Convex Soft Barrier + Urostomy Pouch with Twist Outlet Single Use Kit



- Extended wear barrier + pouch
- Urostomy adapter
- Instructions for use
- Stoma measuring guide

Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
13668	⅝ - 1⅝" (15-40 mm)	
Barrier #16886 Convex Soft, Pouch #122980 Twist outlet / multi-chamber		

Convex Soft Convex Light

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

NOTE: SenSura® Mio Flex and SenSura® Flex **are compatible** with one another

SenSura® Mio Flex Convex Light Barrier + Urostomy Pouch with Twist Outlet Single Use Kit



- Extended wear barrier + pouch
- Urostomy adapter
- Instructions for use
- Stoma measuring guide

Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
13662	⅝ - 1 9⁄16" (15-40 mm)	
Barrier #16483 Convex Light, Pouch #122980 Twist outlet / multi-chamber		

SenSura® Flex Flat Barrier



- Standard wear
- Not made with natural rubber latex

Units: 10
10108† - 5 units

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
10101	⅜ - 1 5⁄16" (10-33 mm)		A4415
10102	⅜ - 1⅞" (10-48 mm)		A4415
10103	⅜ - 2¾" (10-68 mm)		A4415
10108†	⅜ - 3½" (10-88 mm)		A4415

SenSura® Flex Xpro Flat Barrier



- Extended wear
- Not made with natural rubber latex

Units: 10

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
10104	⅜ - 1 5⁄16" (10-33 mm)		A4410
10105	⅜ - 1⅞" (10-48 mm)		A4410
10106	⅜ - 2¾" (10-68 mm)		A4410

Convex Light

SenSura® Flex Convex Light Barrier



- Standard wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex

Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
11301	⅝ - 15⁄16" (15-23 mm)		A4373
11302	⅝ - 1 5⁄16" (15-33 mm)		A4373
11303	⅝ - 1 11⁄16" (15-43 mm)		A4373
11307	⅝ - 2 1⁄16" (15-53 mm)		A4373
Pre-cut			
11312	⅞" (22 mm)		A4373
11313	1" (25 mm)		A4373
11314	1⅛" (28 mm)		A4373
11315	1¼" (31 mm)		A4373

SenSura® Flex Xpro Convex Light Barrier



- Extended wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex

Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
11304	⅝ - 15⁄16" (15-23 mm)		A4408
11305	⅝ - 1 5⁄16" (15-33 mm)		A4408
11306	⅝ - 1 11⁄16" (15-43 mm)		A4408
11308	⅝ - 2 1⁄16" (15-53 mm)		A4408

SenSura® Flex MAXI Drainable Pouch



- EasiClose™ WIDE outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11½" (29½ cm) / 655 mL Units: 20

Code	Code	Color	HCPCS*
Transparent	Opaque		
NA	**11514**		A4425
11512	**11515**		A4425
11517	**11518**		A4425
11½" (29½cm) / 750 mL XXL pouch			
11590	**11591**		A4425

Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

NOTE: SenSura® Mio Flex and SenSura® Flex **are compatible** with one another

SenSura® Flex MAXI Drainable Pouch



- EasiClose™ WIDE outlet
- **Without filter**
- Soft cloth on back
- Transparent
- Not made with natural rubber latex

11½" (29½ cm) / 655 mL — Units: 20

Code	Color	HCPCS*
11612		A5063
11½" (29½cm) / 750 mL XXL pouch		
11690		A5063

SenSura® Flex MAXI Closed Pouch



- **Filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

8½" (21 cm) / 570 mL — Units: 30

Code	Code	Color	HCPCS*
Transparent	Opaque		
10911	**10914**		A4419
10912	**10915**		A4419
10917	**10918**		A4419

Single Use Kits

SenSura® Flex Flat Barrier + MAXI Drainable Pouch Single Use Kit



- Standard wear barrier + pouch
- Instructions for use
- Stoma measuring guide

Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
11011K	⅜ - 1⅞" (10-48 mm)	
Barrier #10102 Flat, Pouch #11512 MAXI / EasiClose™ WIDE outlet / filter		
21001	⅜ - 2¾" (10-68 mm)	
Barrier #10103 Flat, Pouch #11517 MAXI / EasiClose™ WIDE outlet / filter		
21011	⅜ - 3½" (10-88 mm)	
Barrier #10108 Flat, Pouch #11690 XXL / EasiClose™ WIDE outlet / **without filter**		

SenSura® Flex Xpro Flat Barrier + MAXI Drainable Pouch Single Use Kit



- Extended wear barrier + pouch
- Instructions for use
- Stoma measuring guide

Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
21002	⅜ - 1⅞" (10-48 mm)	
Barrier #10105 Flat, Pouch #11512 MAXI / EasiClose™ WIDE outlet / filter		
21003	⅜ - 2¾" (10-68 mm)	
Barrier #10106 Flat, Pouch #11517 MAXI / EasiClose™ WIDE outlet / filter		

SenSura® Flex Xpro Xpro Convex Light Barrier + MAXI Drainable Pouch Single Use Kit



- Extended wear barrier + pouch
- Instructions for use
- Stoma measuring guide

Units: 5

Code	Stoma Size	Color
Cut-to-fit / Transparent		
21004	⅝ - 2 1/16" (15-53 mm)	
Barrier #11308 Convex Light, Pouch #11517 EasiClose™ WIDE outlet / filter		


Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Coloplast

SenSura® Mio Urostomy Pouch Secure with a *twist*





Rosemary
User | SenSura® Mio and Brava®

SenSura® Mio in black

One new color.
One million possibilities.

Having a choice of color for their ostomy pouch reduces the emotional burden of living with an ostomy[1]

Give your patients the power of choice

1. VoxBio study, 2021, Coloplast Market research, Data on file.

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

SenSura® Mio 1-Piece

All the benefits of BodyFit Technology in a flexible low-profile 1-piece system

- 1-piece products have no coupling; the pouch and barrier are welded together as one piece



SenSura® Mio Flat MAXI Drainable Pouch



- Extended wear
- EasiClose™ WIDE outlet
- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile
- Opaque includes inspection window
- Not made with natural rubber latex

11" (28 cm) / 655 mL Units: 10

Code	Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque gray	Opaque black			
10481	**10471**	**18777**	Cut-to-fit	⅜ - 2⅛" (10-55 mm)	A5056
18382	**NA**	**NA**	Cut-to-fit	⅜ - 2¾" (10-70 mm)	A5056
18383	**NA**	**NA**	Cut-to-fit	⅜ - 4" (10-100 mm)	A5056
10482	**10472**	**NA**	Pre-cut	1" (25 mm)	A5056
10488	**10473**	**NA**	Pre-cut	1⅛" (28 mm)	A5056
10483	**10474**	**NA**	Pre-cut	1³⁄₁₆" (30 mm)	A5056
10490	**10475**	**NA**	Pre-cut	1¼" (32 mm)	A5056
10484	**10476**	**NA**	Pre-cut	1⅜" (35 mm)	A5056
10491	**10480**	**NA**	Pre-cut	1½" (38 mm)	A5056
10487	**10479**	**NA**	Pre-cut	1⁹⁄₁₆" (40 mm)	A5056
10485	**10477**	**NA**	Pre-cut	1¾" (45 mm)	A5056
10486	**10478**	**NA**	Pre-cut	2" (50 mm)	A5056

SenSura® Mio Flat MIDI Drainable Pouch



- Extended wear
- EasiClose™ WIDE outlet
- **Full-circle filter**
- Peek-a-boo inspection window allows easy application
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile
- Not made with natural rubber latex

10¼" (26 cm) / 485 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Opaque gray	Opaque black			
10461	**18776**	Cut-to-fit	⅜ - 1¾" (10-45 mm)	A5056

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

SenSura® Mio Flat MAXI Drainable Pouch



- Extended wear
- EasiClose™ WIDE outlet
- **Without filter**
- Water-resistant textile
- Not made with natural rubber latex

Transparent **11" (28 cm) / 655 mL** Units: 10

Code	Cut	Stoma Size	HCPCS*
10489	Cut-to-fit	3/8 - 2 1/8" (10-55 mm)	A4388
10492	Pre-cut	1" (25 mm)	A4388
10493	Pre-cut	1 1/8" (28 mm)	A4388
10494	Pre-cut	1 1/4" (32 mm)	A4388
10495	Pre-cut	1 3/8" (35 mm)	A4388
10496	Pre-cut	1 1/2" (38 mm)	A4388

SenSura® Mio Convex Soft MAXI Drainable Pouch



- Extended wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE outlet
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile
- Opaque includes inspection window
- Not made with natural rubber latex

11" (28 cm) / 655 mL Units: 10

Code	Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque gray	Opaque black			
16715	**16705**	**18990**	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A5057
16716	**16706**	**18991**	Cut-to-fit	3/8 - 2" (10-50 mm)	A5057
NA	**13671**	**NA**	Pre-cut	7/8" (21 mm)	A5057
NA	**13672**	**NA**	Pre-cut	1" (25 mm)	A5057
NA	**13673**	**NA**	Pre-cut	1 1/8" (28 mm)	A5057
NA	**13674**	**NA**	Pre-cut	1 1/4" (31 mm)	A5057
NA	**13675**	**NA**	Pre-cut	1 3/8" (35 mm)	A5057
NA	**13676**	**NA**	Pre-cut	1 1/2" (38 mm)	A5057

Convex Soft

SenSura® Mio Convex Light MAXI Drainable Pouch



- Extended wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE outlet
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile
- Opaque includes inspection window
- Not made with natural rubber latex

11" (28 cm) / 655 mL Units: 10

Code	Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque gray	Opaque black			
16737	**16725**	**18987**	Cut-to-fit	3/8 - 7/8" (10-23 mm)	A5057
16738	**16726**	**18988**	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A5057
16739	**16733**	**18989**	Cut-to-fit	3/8 - 1 11/16" (10-43 mm)	A5057
16740	**16730**	**NA**	Pre-cut	13/16" (21 mm)	A5057
16741	**16731**	**NA**	Pre-cut	1" (25 mm)	A5057
16742	**16732**	**NA**	Pre-cut	1 1/8" (28 mm)	A5057
NA	**16727**	**NA**	Pre-cut	1 1/4" (31 mm)	A5057
NA	**16734**	**NA**	Pre-cut	1 3/8" (35 mm)	A5057
NA	**16735**	**NA**	Pre-cut	1 1/2" (38 mm)	A5057

SenSura® Mio Convex Light MAXI Drainable Pouch



- Extended wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE outlet
- FlexShell is designed for stability and flexibility
- **Without filter**
- Neutral gray, water-resistant textile
- Opaque includes inspection window
- Not made with natural rubber latex

Transparent **11" (28 cm) / 655 mL** Units: 10

Code	Cut	Stoma Size	HCPCS*
16748	Cut-to-fit	3/8 - 1 11/16" (10-43 mm)	A4390

Convex Light

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

SenSura® Mio Convex Deep MAXI Drainable Pouch

- Extended wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE outlet
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile textile with inspection window
- Not made with natural rubber latex



11" (28cm) / 655 mL Units: 10

Code Transparent	Code Opaque gray	Code Opaque black	Cut	Stoma Size	HCPCS*
16765	**16755**	**18984**	Cut-to-fit	3/8 - 7/8" (10-23 mm)	A5057
16766	**16756**	**18985**	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A5057
16767	**16768**	**18986**	Cut-to-fit	3/8 - 1 11/16" (10-43 mm)	A5057
16770	**16760**	**NA**	Pre-cut	13/16" (21 mm)	A5057
16771	**16761**	**NA**	Pre-cut	1" (25 mm)	A5057
16772	**16762**	**NA**	Pre-cut	1 1/8" (28 mm)	A5057
16773	**16757**	**NA**	Pre-cut	1 1/4" (31 mm)	A5057
NA	**16763**	**NA**	Pre-cut	1 3/8" (35 mm)	A5057
NA	**16764**	**NA**	Pre-cut	1 1/2" (38 mm)	A5057

SenSura® Mio Convex Deep MAXI Drainable Pouch

- Extended wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE outlet
- FlexShell is designed for stability and flexibility
- **Without filter**
- Neutral gray, water-resistant textile
- Not made with natural rubber latex



Transparent **11" (28 cm) / 655 mL** Units: 10

Code	Cut	Stoma Size	HCPCS*
16778	Cut-to-fit	3/8 - 1 11/16" (10-43 mm)	A4390

Convex Deep

SenSura® Mio Convex Flip MAXI Drainable Pouch

- Extended wear
- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- Opaque pouches available in neutral gray or black textile with inspection window
- **Full-circle filter**
- Not made with natural rubber latex



11" (28 cm) / 655 mL Units: 10

Code Transparent	Code Opaque gray	Code Opaque black	Cut	Stoma Size	Barrier size**	HCPCS*
18305	**NA**	**NA**	Cut-to-fit	3/8 - 1 5/8" (10-40 mm)	S	A5057
18321	**18312**	**18928**	Cut-to-fit	3/8 - 2" (10-50 mm)	M	A5057
18341	**18330**	**18929**	Cut-to-fit	3/8 - 2 3/8" (10-60 mm)	L	A5057
NA	**18301**	**NA**	Pre-cut	1" (25 mm)	S	A5057
NA	**18302**	**NA**	Pre-cut	1 3/16" (30 mm)	S	A5057
18322	**18313**	**NA**	Pre-cut	1" (25 mm)	M	A5057
18323	**18314**	**NA**	Pre-cut	1 3/16" (30 mm)	M	A5057
18324	**18315**	**NA**	Pre-cut	1 3/8" (35 mm)	M	A5057
18342	**18331**	**NA**	Pre-cut	1 5/8" (40 mm)	L	A5057
18343	**18332**	**NA**	Pre-cut	1 3/4" (45 mm)	L	A5057
18344	**18333**	**NA**	Pre-cut	2" (50 mm)	L	A5057

SenSura® Mio Convex Flip MAXI Drainable Pouch

- Extended wear
- **Without filter**
- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- Not made with natural rubber latex



Transparent **11" (28 cm) / 655 mL** Units: 10


Code	Cut	Stoma Size	Barrier size**	HCPCS*
18320	Cut-to-fit	3/8 - 2" (10-50 mm)	M	A4390
18340	Cut-to-fit	3/8 - 2 3/8" (10–60 mm)	L	A4390

**SenSura Mio Convex Flip is available in multiple barrier sizes to conform to different sizes of outward areas. Start with the smallest barrier size that fits the stoma. For larger outward areas, size up as needed.
- S: Max cut of 1 5/8" (40 mm)
- M: Max cut of 2" (50 mm)
- L: Max cut of 2 3/8" (60 mm)

Convex Flip

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

SenSura® Mio 1-Piece High Output Pouch features



Transparent
- Printed scale for visual indication of output volume

Opaque
- Inspection window
- Option to tuck up outlet for discretion

SenSura® Mio Flat MAXI Drainable Pouch with SOFT Outlet



- Extended wear
- **Full-circle filter designed for ileostomy output**
- For thin, liquid, high-output stool
- Soft outlet can connect to Coloplast® Bed Drainage Bag #14010, or any universal fecal collection bag
- Not made with natural rubber latex

11" (28 cm) / 630 mL — Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque gray			
18665	**18660**	Cut-to-fit	3/8 - 2 3/4" (10-70 mm)	A5056
18666	**18661**	Cut-to-fit	3/8 - 4" (10-100 mm)	A5056

SenSura® Mio Convex Soft MAXI Drainable Pouch with SOFT Outlet



- Extended wear
- Belt tabs for use with 4237, 4247
- **Full-circle filter designed for ileostomy output**
- For thin, liquid, high-output stool
- Soft outlet can connect to Coloplast® Bed Drainage Bag #14010, or any universal fecal collection bag
- Not made with natural rubber latex

11" (28 cm) / 630 mL — Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque gray			
18671	**18670**	Cut-to-fit	3/8 - 2" (10-50 mm)	A5057

SenSura® Mio Convex Light MAXI Drainable Pouch with SOFT Outlet



- Extended wear
- Belt tabs for use with 4237, 4247
- **Full-circle filter designed for ileostomy output**
- For thin, liquid, high-output stool
- Soft outlet can connect to Coloplast® Bed Drainage Bag #14010, or any universal fecal collection bag
- Not made with natural rubber latex

11" (28 cm) / 630 mL — Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque gray			
18676	**18675**	Cut-to-fit	3/8 - 1 11/16" (10-43 mm)	A5057

Convex Soft Convex Light

SenSura® Mio 1-Piece Post Op Pouch with Window



- Neutral gray, water-resistant textile
- Printed scale provides visual indication of output volume
- Easy-to-open and-close post-op window attached to pouch
- Soft outlet can connect to Coloplast® Bed Drainage Bag #14010, or any universal fecal collection bag
- Not made with natural rubber latex

Transparent — Units: 5

Code	Code	Cut	Volume	Stoma Size
Non-sterile	Sterile			
18690	**18680**	Cut-to-fit	680 mL	3/8 - 2 3/4" (10-70 mm)
18691	**18681**	Cut-to-fit	725 mL	3/8 - 4" (10-100 mm)

SenSura® Mio Flat MAXI Closed Pouch



- Standard wear
- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile with inspection window
- Not made with natural rubber latex

8 1/4" (21 cm) / 550 mL — Units: 30

Code	Code	Cut	Stoma Size	HCPCS*
Opaque gray	Opaque black			
10881	**19716**	Cut-to-fit	5/8 - 1 3/4" (15-45 mm)	A4416
10885**	**19719****	Cut-to-fit	1 9/16 - 2 1/8" (40-55 mm)	A4416
10883	**19717**	Pre-cut	1 3/16" (30 mm)	A4416
10884	**19718**	Pre-cut	1 3/8" (35 mm)	A4416
10885**	**19719****	Pre-cut	1 9/16" (40 mm)	A4416

**Because of its larger starter hole, 10885 and 19719 functions as both a pre-cut and cut-to-fit solution.

SenSura® Mio Flat MIDI Closed Pouch



- Standard wear
- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile with inspection window
- Not made with natural rubber latex

7 1/4" (19 cm) / 340 mL — Units: 30

Code	Code	Cut	Stoma Size	HCPCS*
Opaque gray	Opaque black			
10871	**19715**	Cut-to-fit	5/8 - 1 3/4" (15-45 mm)	A4416
10872	**NA**	Pre-cut	1" (25 mm)	A4416
10873	**NA**	Pre-cut	1 3/16" (30 mm)	A4416
10874	**NA**	Pre-cut	1 3/8" (35 mm)	A4416

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

SenSura® Mio Flat MINI Closed Pouch



- Standard wear
- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textilewith inspection window
- Comfortable and discreet stoma cap alternative
- Not made with natural rubber latex

5¼" (13 cm) / 110 mL

Code	Code	Cut	Stoma Size	HCPCS*
Opaque gray	Opaque black			
10870	**19714**	Cut-to-fit	⅝ - 1¾" (15-45 mm)	A4416

SenSura® Mio Convex Soft MAXI Closed Pouch



- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile with inspection window
- Not made with natural rubber latex

8¼" (21 cm) / 550 mL Units: 10

Code	Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque gray	Opaque black			
NA	**16305**	**19701**	Cut-to-fit	⅝ - 1 5/16" (15-33 mm)	A4417
16315	**16306**	**19702**	Cut-to-fit	⅜ - 2" (10-50 mm)	A4417
NA	**16310**	**19703**	Pre-cut	1" (25 mm)	A4417
NA	**16311**	**19704**	Pre-cut	1⅛" (28 mm)	A4417
NA	**16312**	**19705**	Pre-cut	1¼" (31 mm)	A4417
NA	**16313**	**19706**	Pre-cut	1⅜" (35 mm)	A4417

SenSura® Mio Convex Soft MIDI Closed Pouch



- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile with inspection window
- Not made with natural rubber latex

7¼" (19 cm) / 340 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Opaque gray	Opaque black			
16301	**19700**	Cut-to-fit	⅝ - 1 5/16" (15-33 mm)	A4417

Convex Soft

SenSura® Mio Convex Light MAXI Closed Pouch



- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile with inspection window
- Not made with natural rubber latex

8¼" (21 cm) / 550 mL Units: 10

Code	Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque gray	Opaque black			
NA	**16325**	**19573**	Cut-to-fit	⅜ - ⅞" (10-22 mm)	A4417
16335	**NA**	**NA**	Cut-to-fit	⅜ - 1 11/16" (10-43 mm)	A4417
NA	**16326**	**19574**	Cut-to-fit	⅝ - 1 5/16" (15-33 mm)	A4417
NA	**16330**	**19576**	Pre-cut	1" (25 mm)	A4417
NA	**16331**	**19577**	Pre-cut	1⅛" (28 mm)	A4417
NA	**16327**	**19575**	Pre-cut	1¼" (31 mm)	A4417
NA	**16333**	**19579**	Pre-cut	1⅜" (35 mm)	A4417
NA	**16334**	**19580**	Pre-cut	1½" (38 mm)	A4417

SenSura® Mio Convex Light MIDI Closed Pouch



- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile with inspection window
- Not made with natural rubber latex

7¼" (19 cm) / 340 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Opaque gray	Opaque black			
16321	**19571**	Cut-to-fit	⅜ - ⅞" (10-22 mm)	A4417
16322	**19572**	Cut-to-fit	⅝ - 1 5/16" (15-33 mm)	A4417
16323**	**NA**	Cut-to-fit / Pre-cut	1¼ - 1 11/16" (31-43 mm)	A4417

**Because of its larger starter hole, 16323 functions as both a pre-cut and cut-to-fit solution.

Convex Light

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

SenSura® Mio Convex Deep MAXI Closed Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile with inspection window
- Not made with natural rubber latex



8¼" (21 cm) / 550 mL Units: 10

Code Transparent	Code Opaque gray	Code Opaque black	Cut	Stoma Size	HCPCS*
NA	**16345**	**19553**	Cut-to-fit	3/8 - 7/8" (10-22 mm)	A4417
16355	**NA**	**NA**	Cut-to-fit	3/8 - 1 11/16" (10-43 mm)	A4417
NA	**16346**	**19554**	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4417
NA	**16350**	**19556**	Pre-cut	1" (25 mm)	A4417
NA	**16351**	**19557**	Pre-cut	1 1/8" (28 mm)	A4417
NA	**16347**	**19555**	Pre-cut	1 1/4" (31 mm)	A4417
NA	**16353**	**19559**	Pre-cut	1 3/8" (35 mm)	A4417
NA	**16354**	**19560**	Pre-cut	1 1/2" (38 mm)	A4417

SenSura® Mio Convex Deep MIDI Closed Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Water-resistant textile
- Opaque pouches available in neutral gray or black textile with inspection window
- Not made with natural rubber latex



7¼" (19 cm) / 340 mL Units: 10

Code Opaque gray	Code Opaque black	Cut	Stoma Size	HCPCS*
16341	**19551**	Cut-to-fit	3/8 - 7/8" (10-22 mm)	A4417
16342	**19552**	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4417
16343**	**NA**	Cut-to-fit / Pre-cut	1 1/4 - 1 11/16" (31-43 mm)	A4417

**Because of its larger starter hole, 16343 functions as both a pre-cut and cut-to-fit solution.

Convex Deep

SenSura® Mio Convex Flip MAXI Closed Pouch

- Standard wear
- **Full-circle filter**
- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- Opaque pouches available in neutral gray or black textile with inspection window
- Not made with natural rubber latex



8¼" (21 cm) / 550 mL Units: 10

Code Opaque gray	Code Opaque black	Cut	Stoma Size	Barrier size**	HCPCS*
18152	**19532**	Cut-to-fit	3/8 - 1 5/8" (10-40 mm)	S	A4417
18170	**19534**	Cut-to-fit	3/8 - 2" (10-50 mm)	M	A4417
18180	**19538**	Cut-to-fit	3/8 - 2 3/8" (10-60 mm)	L	A4417
18171	**19535**	Pre-cut	1" (25 mm)	M	A4417
18172	**19536**	Pre-cut	1 3/16" (30 mm)	M	A4417
18173	**19537**	Pre-cut	1 3/8" (35 mm)	M	A4417
18181	**19539**	Pre-cut	1 9/16" (40 mm)	L	A4417
18182	**19540**	Pre-cut	1¾" (45 mm)	L	A4417

**SenSura Mio Convex Flip is available in multiple barrier sizes to conform to different sizes of outward areas. Start with the smallest barrier size that fits the stoma. For larger outward areas, size up as needed.
- S: Max cut of 1 5/8" (40 mm)
- M: Max cut of 2" (50 mm)
- L: Max cut of 2 3/8" (60 mm)


Convex Flip

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Urostomy

SenSura® Mio Urostomy Pouch with Twist Outlet



Wide anti-reflux valves support varying ways of managing stents:

- Threading the stents into the pouch via the anti-reflux valves
- Rolling the stents into a loop at the top of the pouch
- Trimming the stents

SenSura® Mio Flat MAXI Urostomy Pouch with Twist Outlet



- Extended wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Twist outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Neutral gray, water-resistant textile on back
- Not made with natural rubber latex

Transparent **10½" (27 cm) / 400 mL** Units: 10 pouches, 5 adapters

Code	Cut	Stoma Size	HCPCS*
105860	Cut-to-fit	⅜ - 1¾" (10-45 mm)	A4428

SenSura® Mio Convex Soft MAXI Urostomy Pouch with Twist Outlet



- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Twist outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Neutral gray, water-resistant textile
- Not made with natural rubber latex

10½" (27 cm) / 400 mL Units: 10 pouches, 5 adapters

Code		Cut	Stoma Size	HCPCS*
Transparent	Opaque gray			
136790	**168050**	Cut-to-fit	⅜ - 1 5/16" (10-33 mm)	A4430
168100	**NA**	Cut-to-fit	⅜ - 2" (10-50 mm)	A4430
137310	**NA**	Pre-cut	⅞" (21 mm)	A4430
137320	**NA**	Pre-cut	1" (25 mm)	A4430
137330	**NA**	Pre-cut	1⅛" (28 mm)	A4430
137340	**NA**	Pre-cut	1¼" (31 mm)	A4430

Convex Soft

SenSura® Mio Convex Soft MIDI Urostomy Pouch with Twist Outlet



- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Twist outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Neutral gray, water-resistant textile
- Not made with natural rubber latex

Opaque gray **9½" (24 cm) / 330 mL** Units: 10 pouches, 5 adapters

Code	Cut	Stoma Size	HCPCS*
168010	Cut-to-fit	⅜ - 1 5/16" (10-33 mm)	A4430

SenSura® Mio Convex Light MAXI Urostomy Pouch with Twist Outlet



- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Twist outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Neutral gray, water-resistant textile
- Not made with natural rubber latex

10½" (27 cm) / 400 mL Units: 10 pouches, 5 adapters

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque gray			
168350	**168250**	Cut-to-fit	⅜ - ⅞" (10-22 mm)	A4430
168360	**168260**	Cut-to-fit	⅝ - 1 5/16" (15-33 mm)	A4430
168370	**NA**	Cut-to-fit	⅜ - 1 11/16" (10-43 mm)	A4430
NA	**168300**	Pre-cut	11/16" (18 mm)	A4430
168400	**168310**	Pre-cut	13/16" (21 mm)	A4430
168410	**168320**	Pre-cut	1" (25 mm)	A4430
168420	**168330**	Pre-cut	1⅛" (28 mm)	A4430
NA	**168270**	Pre-cut	1¼" (31 mm)	A4430

Convex Soft Convex Light

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

SenSura® Mio Convex Light MIDI Urostomy Pouch with Twist Outlet



- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Twist outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Neutral gray, water-resistant textile
- Not made with natural rubber latex

Opaque gray **9½" (24 cm) / 330 mL** Units: 10 pouches, 5 adapters

Code	Cut	Stoma Size	HCPCS*
168210	Cut-to-fit	⅜ - ⅞" (10-22 mm)	A4430
168220	Cut-to-fit	⅝ - 1 5/16" (15-33 mm)	A4430
168230	Pre-cut	1¼" (31 mm)	A4430

SenSura® Mio Convex Deep MAXI Urostomy Pouch with Twist Outlet



- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Twist outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Neutral gray, water-resistant textile
- Not made with natural rubber latex

10½" (27 cm) / 400 mL Units: 10 pouches, 5 adapters

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque gray			
168650	**168550**	Cut-to-fit	⅜ - ⅞" (10-23 mm)	A4430
168660	**168560**	Cut-to-fit	⅝ - 1 5/16" (15-33 mm)	A4430
168670	**NA**	Cut-to-fit	⅜ - 1 11/16" (10-43 mm)	A4430
NA	**168600**	Pre-cut	11/16" (18 mm)	A4430
168700	**168610**	Pre-cut	13/16" (21 mm)	A4430
168710	**168620**	Pre-cut	1" (25 mm)	A4430
168720	**168630**	Pre-cut	1⅛" (28 mm)	A4430
NA	**168570**	Pre-cut	1¼" (31 mm)	A4430

Convex Light Convex Deep

SenSura® Mio Convex Flip MAXI Urostomy Pouch with Twist Outlet



- Extended wear
- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- Twist outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Not made with natural rubber latex

10½" (27 cm) / 400 mL Units: 10 pouches, 5 adapters

Code	Code	Cut	Stoma Size	Barrier size**	HCPCS*
Transparent	Opaque gray				
184050	**184010**	Cut-to-fit	⅜ - 1 5/8" (10-40 mm)	S	A4430
NA	**184020**	Pre-cut	¾" (20 mm)	S	A4430
NA	**184030**	Pre-cut	1" (25 mm)	S	A4430
NA	**184040**	Pre-cut	1 3/16" (30 mm)	S	A4430

**SenSura Mio Convex Flip is available in multiple barrier sizes to conform to different sizes of outward areas. Start with the smallest barrier size that fits the stoma. For larger outward areas, size up as needed.
- S: Max cut of 1 5/8" (40 mm)
- M: Max cut of 2" (50 mm)
- L: Max cut of 2 3/8" (60 mm)

Single Use Kits

SenSura® Mio 1-Piece Flat Single Use Kit



- Extended wear
- Instructions for use
- Stoma measuring guide

Transparent Units: 5

Code	Cut	Stoma Size
21007	Cut-to-fit	⅜ - 2⅛" (10-55 mm)

Pouch #10481 MAXI / EasiClose™ WIDE outlet / filter

SenSura® Mio 1-Piece Convex Soft Single Use Kit



- Extended wear
- Instructions for use
- Stoma measuring guide

Transparent Units: 5

Code	Cut	Stoma Size
13665	Cut-to-fit	⅜ - 2" (10-50 mm)

Pouch #16716 MAXI / EasiClose™ WIDE outlet / filter

Convex Soft Convex Flip

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

SenSura® Mio 1-Piece Convex Light Single Use Kit



- Extended wear
- Instructions for use
- Stoma measuring guide

Transparent — Units: 5

Code	Cut	Stoma Size
13664	Cut-to-fit	⅜ - 1 11/16" (10-43 mm)

Pouch #16739 MAXI / EasiClose™ WIDE outlet / filter

SenSura® Mio 1-Piece Convex Light Urostomy Pouch with Twist Outlet Single Use Kit



- Extended wear + pouch
- Urostomy adapter
- Instructions for use
- Stoma measuring guide

Transparent — Units: 5

Code	Cut	Stoma Size
13663	Cut-to-fit	⅜ - 1 11/16" (10-43 mm)

Pouch #168370 MAXI / Twist outlet / multi-chamber

SenSura® Flat MAXI Drainable Pouch



- Standard wear
- EasiClose™ WIDE outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11½" (29½ cm) / 655 mL — Units: 20

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
15521	**15531**	Cut-to-fit	⅜ - 3" (10-76 mm)	A4424
15522	**15532**	Pre-cut	1" (25 mm)	A4424
15523	**15533**	Pre-cut	1 3/16" (30 mm)	A4424
15524	**15534**	Pre-cut	1⅜" (35 mm)	A4424

SenSura® Flat MAXI Drainable Pouch



- Standard wear
- EasiClose™ WIDE outlet
- **Without filter**
- Soft cloth on back
- Not made with natural rubber latex

Transparent **11½" (29½ cm) / 655 mL** — Units: 20

Code	Cut	Stoma Size	HCPCS*
15821	Cut-to-fit	⅜ - 3" (10-76 mm)	A5061

Convex Light

SenSura® Xpro Flat MAXI Drainable Pouch



- Extended wear
- EasiClose™ WIDE outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11½" (29½ cm) / 655 mL — Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
15981	**NA**	Cut-to-fit	⅜ - 3" (10-76 mm)	A5056
NA	**15991**	Pre-cut	1" (25 mm)	A5056**
NA	**15992**	Pre-cut	1 3/16" (30 mm)	A5056**
NA	**15993**	Pre-cut	1¼" (32 mm)	A5056**
NA	**15994**	Pre-cut	1½" (38 mm)	A5056**

SenSura® Convex Light MAXI Drainable Pouch



- Standard wear
- EasiClose™ WIDE outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11½" (29½ cm) / 655 mL — Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
15606	**15608**	Cut-to-fit	⅝ - 1 11/16" (15-43 mm)	A4389
15623	**15633**	Pre-cut	1" (25 mm)	A4389
15624	**15634**	Pre-cut	1⅛" (28 mm)	A4389

SenSura® Xpro Convex Light MAXI Drainable Pouch



- Extended wear
- EasiClose™ WIDE outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11½" (29½ cm) / 655 mL — Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
15695	**NA**	Cut-to-fit	⅝ - 1 5/16" (15-33 mm)	A5057
15696	**NA**	Cut-to-fit	⅝ - 1 11/16" (15-43 mm)	A5057
NA	**15697**	Pre-cut	⅞" (22 mm)	A5057**
NA	**15698**	Pre-cut	1" (25 mm)	A5057**

Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

**Product has not been through PDAC verification

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

SenSura® Flat MAXI Drainable Pouch



- Standard wear
- Soft outlet provides effective drainage; for thin, liquid, high-output stool; connects to any universal fecal collection bag or Coloplast® Bed Drainage Bag #14010
- **Filter**
- Soft cloth on front and back
- Peek-a-boo inspection window allows easy application
- Not made with natural rubber latex

Opaque **12½" (32 cm) / 630 mL** Units: 30

Code	Cut	Stoma Size	HCPCS*
19004	Cut-to-fit	⅜ - 3" (10-76 mm)	A4424

For urostomy night drainage bag and bag hanger, see page 23.

SenSura® Flat 1-Piece Post-Op Pouch with Window



- Standard wear
- Flat
- Soft outlet can connect to Coloplast® Bed Drainage Bag #14010, or any universal fecal collection bag
- Easy-to-open and - close access window attached to pouch
- Soft cloth on back
- Not made with natural rubber latex

Transparent **13" (33 cm) / 830 mL** Units: 5

Code	Code	Cut	Stoma Size	HCPCS*
Non-sterile	Sterile			
19020	**19010**	Cut-to-fit	⅜ - 3 1/16" (10-78 mm)	A6154
19021	**19011**	Cut-to-fit	⅜ - 4½" (10-115 mm)	A6154

SenSura® 1-Piece Post-Op Pouch



- Standard wear
- Flat
- Soft outlet can connect to Coloplast® Bed Drainage Bag #14010, or any universal fecal collection bag
- Transparent
- Soft cloth on back
- Not made with natural rubber latex

Transparent **12¼" (31 cm) / 680 mL** Units: 5

Code	Cut	Stoma Size	HCPCS*
Non-sterile			
19030	Cut-to-fit	⅜ - 3" (10-70 mm)	A6154

SenSura® 1-Piece Post-Op Pouch with Window



- Standard wear
- Flat
- Soft outlet can connect to Coloplast® Bed Drainage Bag #14010, or any universal fecal collection bag
- Opaque, with inspection window
- Soft cloth on back
- Not made with natural rubber latex

Opaque **12¼" (31 cm) / 630 mL** Units: 5

Code	Cut	Stoma Size	HCPCS*
Sterile			
19001	Cut-to-fit	⅜ - 3" (10-76 mm)	A6154

SenSura® Flat MAXI Closed Pouch



- Standard wear
- **Filter**
- Soft cloth on front and back
- Not made with natural rubber latex

Opaque **8½" (21 cm) / 550 mL** Units: 30

Code	Cut	Stoma Size	HCPCS*
15480	Cut-to-fit	⅜ - 3" (10-76 mm)	A4416
15484	Pre-cut	1⅜" (35 mm)	A4416

Urostomy

SenSura® Flat MAXI Urostomy Pouch



- Standard wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Soft cloth on back
- Not made with natural rubber latex

Transparent **10⅜" (26 cm) / 400 mL** Units: 10 pouches, 5 adapters

Code	Cut	Stoma Size	HCPCS*
11804	Cut-to-fit	⅜ - 3" (10-76 mm)	A5071

SenSura® Flat MIDI Urostomy Pouch



- Standard wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Soft cloth on back
- Not made with natural rubber latex

Transparent **9½" (24 cm) / 330 mL** Units: 10 pouches, 5 adapters

Code	Cut	Stoma Size	HCPCS*
11803	Cut-to-fit	⅜ - 2½" (10-66 mm)	A5071

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

SenSura® Xpro Flat MAXI Urostomy Pouch



- Extended wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

10⅜" (26 cm) / 400 mL Units: 10 pouches, 5 adapters

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
11809	**NA**	Cut-to-fit	⅜ - 3" (10-76 mm)	A4391**
NA	**11805**	Pre-cut	⅞" (22 mm)	A4391**
NA	**11806**	Pre-cut	1" (25 mm)	A4391**
NA	**11807**	Pre-cut	1 3/16" (30 mm)	A4391**
NA	**11808**	Pre-cut	1¼" (32 mm)	A4391**

SenSura® Convex Light MAXI Urostomy Pouch



- Standard wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- Not made with natural rubber latex

Transparent **10⅜" (26 cm) / 400 mL** Units: 10 pouches, 5 adapters

Code	Cut	Stoma Size	HCPCS*
11814	Cut-to-fit	⅝ - 1 5/16" (15-33 mm)	A4392

SenSura® Xpro Convex Light MAXI Urostomy Pouch



- Extended wear
- Belt tabs for use with 4215, 4220
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

10⅜" (26 cm) / 400 mL Units: 10 pouches, 5 adapters

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
11819	**NA**	Cut-to-fit	⅝ - 1 5/16" (15-33 mm)	A4393**
NA	**11816**	Pre-cut	⅞" (22 mm)	A4393**
NA	**11817**	Pre-cut	1" (25 mm)	A4393**
NA	**11818**	Pre-cut	1¼" (32 mm)	A4393**

Indicates convexity

Single Use Kits

SenSura® 1-Piece Flat Single Use Kit



- Standard wear
- Instructions for use
- Stoma measuring guide

Transparent Units: 5

Code	Cut	Stoma Size
11008	Cut-to-fit	⅜ - 3" (10-76 mm)
Pouch #15521 MAXI / EasiClose™ WIDE outlet / filter		

SenSura® Xpro 1-Piece Flat Urostomy Single Use Kit



- Extended wear
- Soft outlet
- Instructions for use
- Stoma measuring guide
- Urostomy adapter

Transparent **10⅜" (26 cm) / 400 mL** Units: 5

Code	Cut	Stoma Size
21005	Cut-to-fit	⅜ - 3" (10-76 mm)
Pouch #11809 MAXI / soft outlet / multi-chamber		

SenSura® Xpro 1-Piece Convex Light Urostomy Single Use Kit



- Extended wear
- Soft outlet
- Instructions for use
- Stoma measuring guide
- Urostomy adapter

Transparent **10⅜" (26 cm) / 400 mL** Units: 5

Code	Cut	Stoma Size
21006	Cut-to-fit	⅝ - 1 5/16" (15-33 mm)
Pouch #11819 MAXI / soft outlet / multi-chamber		

Urostomy Night Bag and Bag Hanger

See Supporting Products page 23 for more information.



Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

**Product has not been through PDAC verification

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

EasiClose™ WIDE outlets offer a 1/2" wider integrated closure for easy, fast emptying on the go.

Assura® Flat MAXI Drainable Pouch



- Standard wear
- EasiClose™ WIDE outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11½" (29 cm) / 530 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
15860	**15870**	Cut-to-fit	⅜ - 3" (10-75 mm)	A4424
NA	**15874**	Pre-cut	1" (25 mm)	A4424
NA	**15875**	Pre-cut	1 3/16" (30 mm)	A4424
NA	**15876**	Pre-cut	1⅜" (35 mm)	A4424

Assura® Flat MAXI Drainable Pouch



- Standard wear
- EasiClose™ WIDE outlet
- **Without filter**
- Soft cloth on front and back
- Not made with natural rubber latex

Opaque **11½" (29 cm) / 530 mL** Units: 10

Code	Cut	Stoma Size	HCPCS*
15871	Cut-to-fit	⅜ - 3" (10-75 mm)	A5061

Assura® Flat MAXI Drainable Pouch



- Extended wear
- EasiClose™ WIDE outlet
- **Without filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11¼" (28½ cm) / 530 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
15869	**15879**	Cut-to-fit	⅜ - 3" (10-75 mm)	A4388

Assura® Convex Light MAXI Drainable Pouch



- Standard wear
- EasiClose™ WIDE outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- Not made with natural rubber latex

Transparent **11½" (29 cm) / 530 mL** Units: 10

Code	Cut	Stoma Size	HCPCS*
14511	Cut-to-fit	⅝ - 1 5/16" (15-33 mm)	A4389
14512	Cut-to-fit	⅝ - 1 11/16" (15-43 mm)	A4389

Assura® Flat MAXI Drainable Pouch



- Standard wear
- EasiClose™ outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11¼" (28½ cm) / 530 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
13860	**13870**	Cut-to-fit	⅜ - 3" (10-75 mm)	A4424
NA	**13874**	Pre-cut	1" (25 mm)	A4424
NA	**13875**	Pre-cut	1 3/16" (30 mm)	A4424
NA	**13876**	Pre-cut	1⅜" (35 mm)	A4424

Assura® Flat MIDI Drainable Pouch



- Standard wear
- EasiClose™ outlet
- **Filter**
- Soft cloth on front and back
- Not made with natural rubber latex

Opaque **10¼" (26 cm) / 390 mL** Units: 10

Code	Cut	Stoma Size	HCPCS*
13840	Cut-to-fit	⅜ - 2½" (10-65 mm)	A4424
13844	Pre-cut	1" (25 mm)	A4424
13845	Pre-cut	1 3/16" (30 mm)	A4424
13846	Pre-cut	1⅜" (35 mm)	A4424

Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Assura® Convex Light MAXI Drainable Pouch



- Standard wear
- EasiClose™ outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11¼" (28½ cm) / 530 mL Units: 10

Code		Cut	Stoma Size	HCPCS*
Transparent	Opaque			
14411	**14414**	Cut-to-fit	⅝ - 1 5⁄16" (15-33 mm)	A4389
14412	**14415**	Cut-to-fit	⅝ - 1 11⁄16" (15-43 mm)	A4389
NA	**14421**	Pre-cut	⅞" (22 mm)	A4389
NA	**14422**	Pre-cut	1" (25 mm)	A4389
NA	**14423**	Pre-cut	1⅛" (28 mm)	A4389
NA	**14424**	Pre-cut	1¼" (31 mm)	A4389

Assura® Convex Light MIDI Drainable Pouch



- Standard wear
- EasiClose™ outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on front and back
- Not made with natural rubber latex

Opaque **10¼" (26 cm) / 390 mL** Units: 10

Code	Cut	Stoma Size	HCPCS*
14404	Cut-to-fit	⅝ - 1 5⁄16" (15-33 mm)	A4389
14405	Cut-to-fit	⅝ - 1 11⁄16" (15-43 mm)	A4389

Assura® Deep Convex MAXI Drainable Pouch



- Standard wear
- EasiClose™ outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11¼" (28½ cm) / 530 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
14103	**14106**	Cut-to-fit	⅝ - 1 11⁄16" (15-43 mm)	A4389
14163	**14173**	Pre-cut	⅞" (22 mm)	A4389
14164	**14174**	Pre-cut	1" (25 mm)	A4389
14165	**14175**	Pre-cut	1⅛" (28 mm)	A4389
14166	**14176**	Pre-cut	1¼" (31 mm)	A4389

Indicates convexity

Assura® Deep Convex MIDI Drainable Pouch



- Standard wear
- EasiClose™ outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on front and back
- Not made with natural rubber latex

Opaque **10¼" (26 cm) / 390 mL** Units: 10

Code	Cut	Stoma Size	HCPCS*
14196	Cut-to-fit	⅝ - 1 11⁄16" (15-43 mm)	A4389

Assura® Convex Light MAXI Drainable Pouch



- Extended wear
- EasiClose™ outlet
- **Without filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11¼" (28½ cm) / 530 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
14413	**NA**	Cut-to-fit	⅝ - 1 5⁄16" (15-33 mm)	A4390
14416	**14419**	Cut-to-fit	⅝ - 1 11⁄16" (15-43 mm)	A4390

Assura® Flat MAXI Drainable Pouch



- Standard wear
- Clamp closure
- **Without filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11¼" (28½ cm) / 530 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
12660	**12670**	Cut-to-fit	⅜ - 2¾" (10-70 mm)	A5061

Assura® Flat MINI Drainable Pouch



- Standard wear
- Clamp closure
- **Without filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

9¾" (24½ cm) / 230 mL Units: 10

Code		Cut	Stoma Size	HCPCS*
Transparent	Opaque			
12630	**12620**	Cut-to-fit	⅜ - 2¼" (10-55 mm)	A5061

Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Assura® Deep Convex MAXI Drainable Pouch



- Standard wear
- Clamp closure
- **Without filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11¼" (28½ cm) / 530 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
13706	**12533**	Cut-to-fit	⅝ - 1 11/16" (15-43 mm)	A4389

Assura® Deep Convex MINI Drainable Pouch



- Standard wear
- Clamp closure
- **Without filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- Not made with natural rubber latex

Transparent **9¾" (24½ cm) / 230 mL** Units: 10

Code	Cut	Stoma Size	HCPCS*
13703	Cut-to-fit	⅝ - 1 11/16" (15-43 mm)	A4389

Assura® Deep Convex MAXI Drainable Pouch



- Extended wear
- Clamp closure
- **Without filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

11¼" (28½ cm) / 530 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
12416	**12419**	Cut-to-fit	⅝ - 1 11/16" (15-43 mm)	A4390

Indicates convexity

Assura® 1-Piece Post-Op Pouch with Window



- Standard wear
- Flat
- Soft outlet can connect to Coloplast® Bed Drainage Bag #14010, or any universal fecal collection bag
- Easy-to-open and -close access window
- Soft cloth on back
- Not made with natural rubber latex

Transparent **13" (33 cm) / 830 mL** Units: 5

Code	Code	Cut	Stoma Size	HCPCS*
Non-sterile	Sterile			
12800	**12808**	Cut-to-fit	⅜ - 2¾" (10-70 mm)	A6154
12802	**12810**	Cut-to-fit	½ - 4" (13-100 mm)	A6154

Assura® 1-Piece Post-Op Pouch

- Standard wear
- Flat
- Soft outlet can connect to Coloplast® Bed Drainage Bag #14010, or any universal fecal collection bag
- Soft cloth on back
- Not made with natural rubber latex

Transparent **12¼" (31 cm) / 690 mL** Units: 5

Code	Code	Cut	Stoma Size	HCPCS*
Non-sterile	Sterile			
12805	**12818**	Cut-to-fit	⅜ - 2¾" (10-70 mm)	A6154

Assura® Flat MAXI Closed Pouch



- Standard wear
- **Filter**
- Soft cloth on front and back
- Not made with natural rubber latex

Opaque **8½" (22 cm) / 530 mL** Units: 30

Code	Cut	Stoma Size	HCPCS*
12170	Cut-to-fit	13/16 - 3" (20-75 mm)	A4416
12174	Pre-cut	1" (25 mm)	A4416
12175	Pre-cut	1 3/16" (30 mm)	A4416
12176	Pre-cut	1⅜" (35 mm)	A4416
12177	Pre-cut	1 9/16" (40 mm)	A4416

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Assura® Flat MIDI Closed Pouch



- Standard wear
- **Filter**
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

7" (18 cm) / 350 mL Units: 30

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
12130	**12140**	Cut-to-fit	13/16 - 2½" (20-65 mm)	A4416
12134	**12144**	Pre-cut	1" (25 mm)	A4416
12135	**12145**	Pre-cut	1 3/16" (30 mm)	A4416
12136	**12146**	Pre-cut	1 3/8" (35 mm)	A4416
NA	**12147**	Pre-cut	1 9/16" (40 mm)	A4416
NA	**12149**	Pre-cut	2" (50 mm)	A4416

Assura® Flat MINI Closed Pouch



- Standard wear
- **Filter**
- Soft cloth on front and back
- Not made with natural rubber latex

Opaque **5" (12 cm) / 170 mL** Units: 30

Code	Cut	Stoma Size	HCPCS*
12110	Cut-to-fit	13/16 - 2½" (20-65 mm)	A4416

Assura® Convex Light MAXI Closed Pouch



- Standard wear
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

8½" (22 cm) / 560 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
14441	**14444**	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4417
14442	**14445**	Cut-to-fit	5/8 - 1 11/16" (15-43 mm)	A4417

Indicates convexity

Assura® Convex Light MIDI Closed Pouch



- Standard wear
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

7" (18 cm) / 350 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
14431	**14434**	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4417
14432	**14435**	Cut-to-fit	5/8 - 1 11/16" (15-43 mm)	A4417
NA	**14451**	Pre-cut	7/8" (22 mm)	A4417
NA	**14452**	Pre-cut	1" (25 mm)	A4417
NA	**14453**	Pre-cut	1 1/8" (28 mm)	A4417
NA	**14454**	Pre-cut	1 1/4" (31 mm)	A4417

Assura® Flat Stoma Cap



- Standard Wear
- **Filter**
- Belt tabs
- Soft cloth on back
- Small pouch size for intimate moments or colostomates who irrigate
- Not made with natural rubber latex

Opaque Units: 30

Code	Cut	Stoma Size	HCPCS*
2501	Cut-to-fit	13/16 - 2¼" (20-55 mm)	A5055

Urostomy

Assura® Flat MAXI Urostomy Pouch



- Standard wear
- Multi-chambers minimize sloshing
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Soft cloth on back
- Not made with natural rubber latex

Transparent **10¾" (27½ cm) / 470 mL** Units: 10 pouches, 5 adapters

Code	Cut	Stoma Size	HCPCS*
14222	Cut-to-fit	3/8 - 2¼" (10-55 mm)	A5071

Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Assura® Convex Light MAXI Urostomy Pouch



- Standard wear
- Multi-chambers minimize sloshing
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- Not made with natural rubber latex

Transparent **10¾" (27½ cm) / 470 mL** Units: 10 pouches, 5 adapters

Code	Cut	Stoma Size	HCPCS*
14717	Cut-to-fit	⅝ - 1 11/16" (15-43 mm)	A4392

Assura® Flat MAXI Urostomy Pouch



- Standard wear
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

10¾" (27 cm) / 465 mL Units: 10 pouches, 5 adapters

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
5585	**5580**	Cut-to-fit	⅜ - 2¼" (10-55 mm)	A5071

Assura® Flat MIDI Urostomy Pouch



- Standard wear
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

7" (18 cm) / 265 mL Units: 10 pouches, 5 adapters

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
5575	**5570**	Cut-to-fit	⅜ - 2¼" (10-55 mm)	A5071

Assura® Flat MAXI Urostomy Pouch



- Extended wear
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Soft cloth on back
- Not made with natural rubber latex

Transparent **9½" (24 cm) / 375 mL** Units: 10 pouches, 5 adapters

Code	Cut	Stoma Size	HCPCS*
12474	Cut-to-fit	⅜ - 2¼" (10-55 mm)	A4391

Indicates convexity

Assura® Flat MIDI Urostomy Pouch



- Extended wear
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

9½" (24 cm) / 265 mL Units: 10 pouches, 5 adapters

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
12472	**12471**	Cut-to-fit	⅜ - 2¼" (10-55 mm)	A4391

Assura® Deep Convex MAXI Urostomy Pouch



- Standard wear
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- Not made with natural rubber latex

Transparent **10¾" (27 cm) / 465 mL** Units: 10 pouches, 5 adapters

Code	Cut	Stoma Size	HCPCS*
12595	Cut-to-fit	⅝ - 1 5/16" (15-33 mm)	A4392
12596	Cut-to-fit	⅝ - 1 11/16" (15-43 mm)	A4392
12991	Pre-cut	⅝" (15 mm)	A4392
12992	Pre-cut	¾" (18 mm)	A4392
12993	Pre-cut	⅞" (22 mm)	A4392
12994	Pre-cut	1" (25 mm)	A4392
12995	Pre-cut	1⅛" (28 mm)	A4392
12996	Pre-cut	1¼" (31 mm)	A4392
12997	Pre-cut	1⅜" (35 mm)	A4392

Assura® Deep Convex MAXI Urostomy Pouch



- Extended wear
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- Not made with natural rubber latex

Transparent **10¾" (27 cm) / 465 mL** Units: 10 pouches, 5 adapters

Code	Cut	Stoma Size	HCPCS*
12432	Cut-to-fit	⅝ - 1 5/16" (15-33 mm)	A4393
12433	Cut-to-fit	⅝ - 1 11/16" (15-43 mm)	A4393

Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Single Use Kits




Assura® 1-Piece Flat Single Use Kit

- Standard wear
- Instructions for use
- Stoma measuring guide

Transparent — Units: 5

Code	Cut	Stoma Size
11005	Cut-to-fit	⅜ - 3" (10-76 mm)

Pouch #15860 MAXI / EasiClose™ WIDE outlet / filter



Coloplast® Fistula and Wound Management System

Each system includes:

- 1 pouch and tracing guide
- 1 drain port (extras sold separately on next page)
- 1 lid (extras sold separately on next page)
- 1 clear drape
- 1 clamp (extras sold separately on next page)
- Not made with natural rubber latex

Each box includes one pump.



14050 Mini



14060 Midi



14070 Maxi

Units: 3
Units for 14050: 6

Code	Cut	Stoma Size	HCPCS*
14050	Mini	4⅛ x 6¼" (104-159 mm)	A6154
14060	Midi	6⅛ x 9" (156-228 mm)	A6154
14070	Maxi	8⅛ x 11¾" (208-297 mm)	A6154

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Fistula / Wound Pouch Accessories

Coloplast® Fistula Drain Port



For use with Fistula and Wound Management Systems. Drain port can be easily applied to the lid to allow for catheter insertion.

- For catheters ranging from 18FR - 22FR, the drain port is ready for use
- For catheters ranging from 24FR - 30FR, pull off the top of the drain port -DO NOT CUT
- Not made with natural rubber latex

Code	Size	Units
14015	One size	12

Coloplast® Bed Drainage Bag



Wide-lumen bag tubing specifically designed for fecal output.

- Holds up to 2 liters (2000 mL)
- Can be cut to find right length for bed drainage
- Compatible with fistula and wound management systems or high-output pouches
- Bag hanger (#5070) not included- order separately
- Not made with natural rubber latex

Units: 6

Code	Size	HCPCS*
14010	2 liters (2000 mL)	A4357

For matching bag hanger, see page 23.

Coloplast® Transparent Fistula Lid



Replacement lids match the transparent lid included with each fistula pouch.

- Not made with natural rubber latex

Code	Size	Units
14021	Mini (matches 14050)	10
14031	Midi (matches 14060)	10
14041	Maxi (matches 14070)	10

Coloplast® Pouch Clamp



Simple closure clamp is easy to manipulate, even for those with low dexterity.

- Not made with natural rubber latex

Units: 20

Code	Size	HCPCS*
9500	One size	A4363

SenSura® Mio Baby & Kids


For the little fighters

Babies and kids present different pouching challenges than adults - so we've developed specific solutions for different stages in a child's life.





Stepping up over time, as kids grow



Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Ask us about our pediatric education!




Pediatric FAQ for parents & caregivers
Tipsheet with answers to frequently asked questions, like "what do I do if friends and family - or other kids - have questions?" and "what do I do if my child pulls the pouch off?" Available in English and Spanish.
Caregivers can request this from our Coloplast Care program.



My child has an ostomy: resource sheet for caregivers
A fill-in-the-blank sheet to help parents communicate about their child's needs to daycare providers, school nurses, and other caregivers.
Caregivers can request this from our Coloplast Care program.



Life after your child's ostomy
A booklet covering an overview of anatomy, how to change a pouch, pouching system options, where to get supplies, dietary notes, and other information. Available in English and Spanish.
Caregivers can request this from our Coloplast Care program.



When I met Tipster: A coloring book for children with ostomies
A simple coloring book to help kids engage with their ostomy surgery and "new normal" - or learn more about a sibling's ostomy.
Available to healthcare facilities via Coloplast sales representatives, or caregivers can request from our Coloplast Care program.




Jake Gets a Stoma
A children's book written by Jake's mom, who wanted a better way to explain the surgery and stoma management to her son and other kids going through the same thing.
Available to healthcare facilities via Coloplast sales representatives, or caregivers can request from our Coloplast Care program.



SenSura® Mio Bunny
Coloplast offers a complimentary training bunny for product education use with both SenSura® Mio Baby & Kids. *Available to healthcare facilities via Coloplast sales representatives.*

Limitations apply. The SenSura® Mio Bunny is provided to healthcare facilities for families for teaching purposes only and should only be used with Coloplast products.



Teen life with an ostomy
Navigating school, what to wear, and hanging out with friends aren't covered by the average adult-oriented materials - that's why we're so excited to feature tips from Paige, a teen who's been through it all!
Caregivers can request this from our Coloplast Care program.

Information from Coloplast Care is for educational purposes only. It is not intended to substitute for professional medical advice and should not be interpreted to contain treatment recommendations. You should rely on the healthcare professional who knows your individual history for personal medical advice and diagnosis.

SenSura® Mio Baby

2-piece pouching system designed to tackle the challenges of preemies and babies from 30 weeks of gestation to 6 months of age.

SenSura® Mio Baby Flex 2-Piece Barrier



- Oval barrier accommodates one or multiple stomas
- No starter hole – custom cut
- Gentle adhesive sticks securely
- Removes gently with Brava® Baby Adhesive Remover
- Soft flange enables the pouch to be removed and replaced
- Sterilized surface, single-packed
- Not made with natural rubber latex

Cut-to-fit — Units: 5

Code	Stoma Size	Color	HCPCS*
18705	0 - 19/16" (0-40 mm)		A4414

SenSura® Mio Baby Flex 2-Piece Pouch



- Adhesive coupling is easy to apply and remove
- Access window to refeed or expel gas
- Soft outlet fits to enteral feeding systems
- Sterile, single-packed
- Not made with natural rubber latex

Transparent, No filter — Units: 10

Code	Color	HCPCS*
18700		A5063

Brava® Baby Adhesive Remover



- Oil-based adhesive remover designed specifically for babies
- Liquid sachet (no wipe) can be used with cotton swabs or gauze
- Gently removes adhesive and adhesive residues from SenSura® Mio Baby Flex 2-Piece Barrier
- Not recommended for use with SenSura® Mio Kids or adult products
- Not made with natural rubber latex

Units: 30

Code	Size	HCPCS*
12014	3 mL sachet	A4455

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

SenSura® Mio Kids

1-piece and 2-piece pouching systems designed with BodyFit Technology to fit kids' active movements and rounded bellies. Designed with options to step up as kids grow, the range covers kids 6 months to 4 years of age (larger options for older kids are found in the SenSura Mio adult line).

2-Piece

SenSura® Mio Kids Flex 2-Piece Barrier



- Star-shaped, elastic barrier can be rotated for flexibility in placement
- No starter hole – custom cut
- Not made with natural rubber latex

Cut-to-fit — Units: 5

Code	Stoma Size	Color	HCPCS*
18735	0 - 1 1/16" (0-27 mm)		A4414
18747	0 - 1 3/8" (0-35 mm)		A4414

SenSura® Mio Kids Flex 2-Piece MINI Drainable Pouch



- EasiClose™ outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Small size enables kids to step up in pouch size over time
- Opaque gray includes inspection window
- Not made with natural rubber latex

6 1/8" folded — Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque gray		
18721	**18720**		A4425

SenSura® Mio Kids Flex 2-Piece MIDI Drainable Pouch



- EasiClose™ outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque gray includes inspection window
- Not made with natural rubber latex

6¾" folded — Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque gray		
18723	**18722**		A4425
18725	**18724**		A4425

SenSura® Mio Kids Flex 2-Piece MIDI Urostomy Pouch



- Neutral gray, water-resistant textile
- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for kids in a wheelchair
- Connects to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag with enclosed adapters
- Not made with natural rubber latex

Opaque gray **6¾" folded** — Units: 10 pouches, 5 adapters

Code	Color	HCPCS*
18730		A5073
18731		A5073

1-Piece

SenSura® Mio Kids 1-Piece MIDI Drainable Pouch



- Star-shaped, elastic barrier with no starter hole – custom cut
- EasiClose™ outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque gray includes inspection window
- Not made with natural rubber latex

6¾" folded | Cut-to-fit — Units: 10

Code	Code	Stoma size	HCPCS*
Transparent	Opaque gray		
18711	**18710**	0 – 1 3/8" (0-35 mm)	A4424

SenSura® Mio Kids 1-Piece MIDI Urostomy Pouch



- Star-shaped, elastic barrier with starter hole
- Neutral gray, water-resistant textile
- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for kids in a wheelchair
- With enclosed adapters, connects to Coloplast® Urostomy Night Bag and any universal bedside drainage bag
- Not made with natural rubber latex

Cut-to-fit — Units: 10 pouches, 5 adapters

Code	Code	Stoma size	HCPCS*
Transparent	Opaque gray		
18712	**18713**	3/8 - 1 3/8" (10-35 mm)	A5071

Urostomy Night Bag and Bag Hanger

See Supporting Products on page 23 for more information.





*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Assura® AC (Adhesive Coupling) Pediatric

A simple, soft, secure system for neonates. The pouching system is easy to clean, and simple to handle with the EasiClose™ standard integrated outlet.

Assura® AC Pediatric Flat Barrier



- Standard wear
- No starter hole
- Not made with natural rubber latex

Cut-to-fit — Units: 10

Code	Stoma Size	Color	HCPCS*
14307	0 - ⅝" (0-17 mm)		A4414
14308	0 - 1" (0-27 mm)		A4414

Assura® AC Pediatric Drainable Pouch



- EasiClose™ outlet
- **Filter**
- Soft cloth on back
- Not made with natural rubber latex

Units: 10

Code	Code	Length	Color	HCPCS*
Transparent	Opaque			
14681	**14691**	5½" (14 cm)		A4425
14682	**14692**	6½" (16 cm)		A4425

Assura® 2-Piece Pediatric Pouching Systems

Specifically developed adhesive and a smaller barrier size for pediatrics, while a mix-and-match "green" flange size offers versatility. The SenSura® and Assura® pouches and barriers may be used interchangeably with "green"-sized adult SenSura® Click or Assura® products for a broad array of pouching options.

Assura® Pediatric Standard Wear Barrier



- Flat
- Not made with natural rubber latex

Cut-to-fit — Units: 5

Code		Color	HCPCS*
2182	⅜ - 1⅜" (10-35 mm)		A4414

Assura® Pediatric Drainable Pouch



- Clamp closure
- Soft cloth on back
- Not made with natural rubber latex

Transparent **8½" (21 cm) / 150 mL** — Units: 10

Code	Color	HCPCS*
2156		A4426

Assura® Pediatric Closed Pouch



- Soft cloth on front and back
- Not made with natural rubber latex

Opaque **5¾" (14½ cm)** — Units: 30

Code	Color	HCPCS*
2161		A4423

Assura® Stoma Cap



- **With filter**
- Small size is optimal for swimming at the pool
- Not made with natural rubber latex

Opaque — Units: 30

Code	Color	HCPCS*
2801		A5055

Assura® Pediatric Urostomy Pouch



- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter - see page 23 for options
- Soft cloth on back
- Not made with natural rubber latex

Transparent **8½" (21 cm) / 150 mL** — Units: 10 pouches, 5 adapters

Code	Color	HCPCS*
1789		A4433

Assura® Urostomy Micro-Pouch



- Miniature size prevents twisting
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter - see page 23 for options
- Soft cloth on back
- Not made with natural rubber latex

Transparent **6" (15 cm)** — Units: 10

Code	Color	HCPCS*
14204		A4433

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Assura® Urostomy Minicap



- Ultra-absorbent material inside
- Small size is optimal for swimming at the pool
- Not made with natural rubber latex

Opaque **100 mL** — Units: 30

Code	Color	HCPCS*
2807		A5055

ColoKids™ 1-Piece Pediatric Pouching Systems

Utilizes Assura® technology and features small sizes for pediatric-friendly pouching options.

ColoKids™ Flat 1-Piece Drainable Pouch



- Standard wear
- Clamp closure
- Soft cloth on back
- Not made with natural rubber latex

Transparent — Units: 10

Code	Cut	Stoma Size	HCPCS*
2115	Cut-to-fit	⅜ - 1⅜" (10-35 mm)	A5061

ColoKids™ Flat 1-Piece Closed Pouch



- Standard wear
- Soft cloth on back
- Not made with natural rubber latex

Transparent — Units: 30

Code	Cut	Stoma Size	HCPCS*
2125	Cut-to-fit	⅜ - 1⅜" (10-35 mm)	A4416

ColoKids™ Flat 1-Piece Urostomy Pouch



- Standard wear
- Soft outlet can be connected to Coloplast® Urostomy Night Bag #21365 and any universal bedside drainage bag or leg bag using included adapter - see page 23 for options
- Soft cloth on back
- Not made with natural rubber latex

Transparent — Units: 10 pouches, 5 adapters

Code	Cut	Stoma Size	HCPCS*
8009	Cut-to-fit	⅝ - 1⅜" (15-35 mm)	A5071

Urostomy Night Bag and Bag Hanger

See Supporting Products on page 23 for more information.





Brava®

Using Brava® Supporting Products with Pediatrics

Based on updated documentation guidelines, real-life use, and existing biological evaluation reports, Coloplast reassessed the risk to pediatrics with the use of supporting products. Updated guidelines are provided below.

Product	Caution
Brava® Protective Seal	No caution
Brava® Protective Seal Convex	
Brava® Elastic Barrier Strips	
Brava® Protective Sheet	
Brava® Lubricating Deodorant	
Brava® Skin Cleanser Wipes	
Brava® Moldable Ring	*This product is non-sterile and should be used with caution on premature babies or infants receiving incubator care as it may lead to increased risks of infection.*
Brava® Strip Paste	
Brava® Paste	
Brava® Powder	
Brava® Adhesive Remover Wipes	*Not recommended to be used on infants and neonates below 7 lbs 11 oz (3.5 kg) because of their immature skin.*
Brava® Skin Barrier Wipes	
Brava® Adhesive Remover Spray	*Not recommended to be used on infants and neonates below 7 lbs 11 oz (3.5 kg) because of their immature skin. Avoid direct inhalation as it may give rise to temporary irritation of the airways.*
Brava® Skin Barrier Spray	

Any use outside of our indication is considered off-label use and is up to the clinician's own discretion. Clinicians must weigh the risks and benefits of available products to make an informed and appropriate decision for each patient.

*For important information regarding the use of HCPCS codes listed and code definitions see page 124

SenSura® Mio carries patients through every stage of life

SenSura® Mio Baby & Kids





SenSura® Mio High Output and Post-Op





SenSura® Mio









Because *every fit* matters



Making life easier_



Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Coloplast | Care

Delivering Personalized Care Solutions.

Simply. Reliably. For You.

About the team



Coloplast Care Advisors assist patients with their intimate healthcare needs, by providing personalized lifestyle and product advice and ongoing resources. We will assist in coordinating communications with a supplier to place an order when appropriate, or refer back to clinicians if there are concerns.

Tell your patients we will be calling from:
1-877-858-2656 Hablamos español
www.ostomy.coloplastcare.us

Ostomy patients enrolled experienced:

~30%↓ Reduction in hospital readmission rates†1

~45%↓ Reduction in unnecessary ER visits‡1

~25%↑ Increase in HCAHPS satisfaction scores†2

Each person's situation is unique and risks, outcomes, experience, and results may vary.

Ostomy patients that engage with our Coloplast Care Advisors do better over time and had statistically fewer peristomal skin conditions[3]

†Observed values within the first month of discharge.
‡Observed values after the first month of discharge.

1. Rojanasarot, S., The Impact of a Readily Available, Post-Discharge Support Program for Ostomy Surgery Patients on Preventable Health Care Utilization, Jan/Feb 2018 JWOCN.
2. Helder, D., O'Day, M., Measuring System Impact by Improving the Ostomy Patient Experience from a Product Standardization Project in a 13-Hospital Integrated Delivery Network in the Midwest. J WOCN, 2017, 44(35): S53
3. Burton, M, et. al, The positive Impact of Coloplast Care, a Patient-support program: results from Patient reported outcomes, Poster Presentation, NorthEast Regional Affiliate of the WOCN Society, Oct 2023.

Information from Coloplast Care is for educational purposes only. It is not intended to substitute for professional medical advice and should not be interpreted to contain treatment recommendations. You should rely on the healthcare professional who knows your individual history for personal medical advice and diagnosis.




Coloplast | Care

Coloplast Care *Connect*

***Connect* can help you:**

- Streamline ordering and enrollment
- View a dashboard of enrolled patients within your facility
- Share workflows and finish documentation before the patient transitions



Accessible on desktop and mobile devices

www.coloplastcareconnect.us

Scan the QR code with your smartphone camera and select **"New Enrollment"** to get started!

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services



Trust our Clinician Support Line to support your discharge process for ostomy patients

Call us at 1-855-863-3912

Our Care Advising teams participate in industry education from WEB WOC® Programs, an accredited and internationally recognized WOC education provider.

WEB WOC is a registered trademark of WEB WOC nursing education program. Used with permission.

Clinician Support Line

How does the discharge process for ostomy patients work?

Patient is enrolled




Coloplast Care Welcome Kit* ships to patient's home



Advisor calls patient to educate and confirm insurance for ordering supplies



Patient is connected with a supplier for fulfillment



Supplier ships and bills insurance



*Limitations apply

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Coloplast Care Welcome Kits are tailored to include Living Well books specific to each type of ostomy, along with appropriate items sent per ostomy type.

Ostomy Self-Assessment Tools
- BodyCheck
- Ostomy Check
- Troubleshooter
- Apply-Remove-Check (ARC)

Diaries:
- Pouch change diary
- Exercise diary
- Hernia prevention action plan
- Ostomy Wellness Output Diary

Educational Guides
- Life with your ostomy (ileostomy, colostomy, urostomy)
- Life after your child's ostomy
- Teen life with an ostomy

Product Pictorial Guide - Also available in Spanish
NOTE: These are abbreviated instructions designed to visually show how to use the product. Always refer to the product labeling (including the package insert) before using.
- SenSura® Mio 1-Piece Pouching System
- SenSura® Mio Click 2-Piece Pouching System
- SenSura® Mio Flex 2-Piece Pouching System
- SenSura® 1-Piece Pouching System
- SenSura® Click 2-Piece Pouching System
- SenSura® Flex 2-Piece Pouching System
- SenSura® Mio 1-Piece Urostomy Pouching System
- SenSura® Mio Click 2-Piece Urostomy Pouching System
- SenSura® Mio Flex 2-Piece Urostomy Pouch System
- SenSura® Mio Baby 2-Piece Flex Barrier, Pouch, and Brava® Baby Adhesive Remover

Spanish Educational Tools
Educational Guides
- Life with your ostomy (ileostomy, colostomy, urostomy)
- Life after your child's ostomy

Also please note, that our Coloplast Care website can be translated in Chrome browser, **www.ostomy.coloplastcare.us**

Other languages also available in digital format only: Arabic, Chinese, and Russian. These can be found at the link below.

For additional educational resources, visit
https://coloplast.to/patienteducation

Coloplast | Professional

Supporting your care of every patient

You play an important role in the lives of your patients, which is why we want to support your professional development and day-to-day work for your patient's benefit and for yours.

www.coloplastprofessional.us



Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

SenSura® Mio 2-piece pouching system



Color
(Barrier and pouch must use matching color products)





A colored circle at the bottom of the barrier matches the box color
(Barrier product number is shown on the box label)





Product number
(Can be found by lifting the fabric at the top of the pouch)

SenSura® Mio 1-piece pouching system







Product number
(Can be found by lifting the fabric at the top of the pouch)



Convexity type

Identifying convexity

From the box

Convexity type - *Check product guide for availability*

Soft	Light	Deep	Flip
Most compressible		Most stable (least compressible)	

Identifying convexity

From the item

- Convex barriers have a dome-like structure
- Convex Flip barriers have a star shape



Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Prior to use, refer to the product 'Instructions for use' for intended use and relevant safety information.

Cut the barrier to fit the exact size and shape of the stoma





Make sure the skin around the stoma is clean and dry



- Use water and a clean cloth
- Ensure all residue is removed
- Do not use soaps, moisturizers, or oil-based products as these may affect adhesion
- Consult your Wound Ostomy Continence (WOC) nurse for questions about supporting products

Apply gentle hand pressure to the entire barrier to initiate adhesion




- Run a finger around the stoma, then smooth outwards to the edges
- Ensure the entire barrier has been massaged into the skin

For product support call the Coloplast Care Clinician Support Line at 1-855-863-3912

Prior to use, refer to the product 'Instructions for use' for intended use and relevant safety information.

Close drainable pouches prior to use



- For a drainable fecal pouch, disconnect the VELCRO® brand dots used to empty, then fold the outlet towards the head of the patient three times. Fold the two turqoise tabs inwards to secure the VELCRO® brand fastener.

When using a 2-piece system, grab both a barrier and a pouch with a matching color code

- For **SenSura® Mio Flex**, align the edge of the adhesive ring with the turquoise ring on the outer rim of the barrier and apply gentle pressure around the barrier coupling to secure.


SenSura® Mio Flex

- For **SenSura® Mio Click**, align the coupling rings and, keeping your hand flat against the pouch, gently press the pouch to the barrier – you should feel it 'click' into place. Squeeze the lock to secure the pouch.


SenSura® Mio Click

For product support call the Coloplast Care Clinician Support Line at 1-855-863-3912

VELCRO® is a registered trademark of Velcro BVBA. Used with permission.

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Medicare Utilization - Ostomy

Usual maximum limits for some products from the Medicare Local Coverage Determinations (LCD)[1] are shown below. The amount of ostomy supplies varies by patient based on the type of ostomy, its location, its construction, and the condition of the peristomal skin. The table below shows the maximum number of items that are usually medically necessary for patients with ostomies. It may be more or less depending on the factors that affect the frequency of barrier and pouch change. The medical necessity for a greater quantity of supplies must be clearly documented in the patient's medical record (physician's signed and dated order).

HCPCS	Short Description	Usual Maximum[1]
A4357	Bedside drainage bag	2 ea/month
A4361	Ostomy faceplate, each	3 ea/6 months
A4362	Skin barrier, solid, 4x4 or equivalent; each**	20 ea/month
A4363	Ostomy clamp replacement	none*
A4367	Ostomy belt, each	1 ea/month
A4368	Ostomy filter	none*
A4369	Ostomy skin barrier; liquid, per ounce	2 oz./month
A4371	Ostomy skin barrier powder, per ounce	10 oz/6 months
A4372	Ostomy skin barrier, 4 x 4, standard wear, CONVEX	none*
A4373	Ostomy skin barrier, with flange, CONVEX	none*
A4385	Ostomy skin barrier, 4 x 4, extended wear, FLAT	none*
A4388	One-piece ostomy drainable pouch and barrier, extended wear, FLAT	none*
A4389	One-piece ostomy drainable pouch and barrier, CONVEX	none*
A4390	One-piece ostomy drainable pouch, extended wear barrier, CONVEX	none*
A4391	One-piece urinary pouch, extended wear barrier, FLAT	none*
A4392	One-piece urinary pouch, standard wear barrier, FLAT	none*
A4394	Ostomy deodorant per ounce	none*
A4396	Ostomy belt with peristomal hernia support	none*
A4398	Ostomy irrigation supply; bag, each	2 ea/6 months
A4399	Ostomy irrigation supply; cone/catheter, with or without brush	2 ea/6 months
A4404	Ostomy ring, each	10 ea/month
A4405	Ostomy skin barrier, non-pectin based, paste, per ounce	4 oz/month
A4406	Ostomy skin barrier, pectin-based, per ounce	4 oz/month
A4407	Two-piece ostomy skin barrier, with flange, extended wear, CONVEX, (4 x 4 inches or smaller)	none*
A4408	Two-piece ostomy skin barrier, with flange extended wear, CONVEX, (larger than 4 x 4 inches)	none*
A4409	Two-piece ostomy skin barrier, with flange, extended wear, FLAT, (4 x 4 inches or smaller)	none*
A4410	Two-piece ostomy skin barrier, with flange, extended wear, FLAT, (larger than 4 x 4 inches)	none*
A4411	Ostomy skin barrier, solid 4x4 or equivalent, extended wear, with built-in convexity, each	none*
A4412	Two-piece high output ostomy drainable pouch without filter	none*

* No published Medicare usual maximum. Monthly quantity subject to medical necessity documentation.
**Elastic Barrier Strips - Curved or Straight: 40 strips/month

Medicare utilization – ostomy

HCPCS	Short Description	Usual Maximum[1]
A4413	Two-piece high output ostomy drainable pouch with filter	none*
A4414	Two-piece ostomy skin barrier with flange, FLAT (4x4 inches or smaller)	20 ea/month
A4415	Two-piece ostomy skin barrier, with flange, FLAT, greater than 4 x 4 inches	20 ea/month
A4416	One-piece ostomy closed pouch and barrier with filter, FLAT	60 ea/month
A4417	One-piece ostomy closed pouch and barrier with filter, CONVEX	60 ea/month
A4418	One-piece ostomy closed pouch with filter and barrier	60 ea/month
A4419	Two-piece ostomy closed pouch with filter for non-locking system	60 ea/month
A4423	Two-piece ostomy closed pouch with filter for locking system	60 ea/month
A4424	One-piece ostomy drainable pouch and barrier with filter, FLAT	20 ea/month
A4425	Two-piece ostomy drainable pouch with filter for non-locking system	20 ea/month
A4426	Two-piece ostomy drainable pouch without filter for locking system	20 ea/month
A4427	Two-piece ostomy drainable pouch with filter for locking system	20 ea/month
A4433	Two-piece ostomy urinary pouch for locking system	20 ea/month
A4436	Irrigation supply, sleeve, disposable, per month	1 month supply
A4437	Irrigation supply, sleeve, disposable, per month	1 month supply
A4455	Adhesive remover or solvent, per ounce	16 oz/6 months
A4456	Adhesive remover, wipes	none*
A5055	Stoma cap	31 ea/month
A5056	One-piece drainable ostomy pouch, extended wear with filter	40 ea/month
A5057	One-piece drainable ostomy pouch extended wear with filter, CONVEX	40 ea/month
A5061	One-piece ostomy drainable pouch and barrier, FLAT	20 ea/month
A5063	Two-piece ostomy drainable pouch for barrier with flange	20 ea/month
A5071	One-piece ostomy urinary pouch	20 ea/month
A5073	Two-piece ostomy urinary pouch	20 ea/month
A5120	Skin barrier, wipes or swabs, each	150 ea/6 months
A5121	Solid skin barrier, 6x6 or equivalent, each**	20 ea/month
A5122	Solid skin barrier, 8x8 or equivalent, each**	20 ea/month
A6154	Wound pouch, each	3/week[2]

* No published Medicare usual maximum. Monthly quantity subject to medical necessity documentation.
**Elastic Barrier Strips - Y-shape or XL: 60 strips/month

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

Reimbursement Disclaimer: Coloplast Corp. provides this information for convenience and your general reference only. It does not constitute legal advice or a recommendation regarding clinical practice. Reimbursement, coverage and payment policies can vary from one insurer and region to another and is subject to change without notice. The provider has the responsibility to determine medical necessity and to submit appropriate codes and charges for care provided. Coloplast does not guarantee coverage or payment of products and Coloplast makes no guarantee that the use of this information will prevent differences of opinion or disputes with Medicare, insurers, or other payers as to the correct form of billing or the amount that will be paid. This information is provided for your general information only and is not intended to replace any advice you receive from your own internal or external insurance coverage consultants, reimbursement specialists or legal counsel.

HCPCS Code Reference	
A4357	Bedside drainage bag, day or night, with or without anti-reflux device, with or without tube, each
A4361	Ostomy faceplate, each
A4362	Solid skin barrier; solid, 4x4 or equivalent; each
A4363	Ostomy clamp, any type, replacement only, each
A4366	Ostomy vent, any type, each
A4367	Ostomy belt, each
A4368	Ostomy filter, any type, each
A4369	Ostomy skin barrier liquid (spray, brush, etc.), per ounce
A4371	Ostomy skin barrier powder, per ounce
A4372	Ostomy skin barrier, solid 4x4 or equivalent, standard wear, with built-in convexity, each
A4373	Ostomy skin barrier, with flange (solid, flexible or accordion), with built-in convexity, any size, each
A4385	Ostomy skin barrier, solid 4 x 4 or equivalent, extended wear, without built-in convexity, each
A4388	Ostomy pouch, drainable, with extended wear barrier attached, (1 piece), each
A4389	Ostomy pouch, drainable, with barrier attached, with built-in convexity (1-piece), each
A4390	Ostomy pouch, drainable, with extended wear barrier attached, with built-in convexity (1-piece), each
A4391	Ostomy pouch, urinary, with extended wear barrier attached (1-piece), each
A4392	Ostomy pouch, urinary with standard wear barrier attached, with built-in convexity (1-piece), each
A4393	Ostomy pouch, urinary, with extended wear barrier attached, with built-in convexity (1 pc), each
A4394	Ostomy deodorant, with or without lubricant, for use in ostomy pouch, per fluid ounce
A4396	Ostomy belt with peristomal hernia support
A4398	Ostomy irrigation supply; bag, each
A4399	Ostomy Irrigation supply; cone /catheter, with or without brush
A4404	Ostomy ring, each
A4405	Ostomy skin barrier, non-pectin based, paste, per ounce
A4406	Ostomy skin barrier, pectin-based, paste, per ounce
A4407	Ostomy skin barrier, with flange (solid, flexible, or accordion), extended wear, with built-in convexity, 4x4 inches or smaller, each
A4408	Ostomy skin barrier, with flange (solid, flexible or accordion), extended wear, with built-in convexity, larger than 4x4 inches, each

HCPCS code reference

HCPCS Code Reference	
A4409	Ostomy skin barrier, with flange (solid, flexible or accordion), extended wear, without built in convexity, 4x4 inches or smaller, each
A4410	Ostomy skin barrier with flange (solid, flexible or accordion), extended wear, without built-in convexity, larger than 4x4 inches, each
A4411	Ostomy skin barrier, solid 4x4 or equivalent, extended wear, with built-in convexity, each
A4412	Ostomy pouch, drainable, high output, for use on a barrier with flange (2-piece system), without filter, each
A4413	Ostomy pouch, drainable, high output, for use on a barrier with flange (2 piece system) with filter, each
A4414	Ostomy skin barrier, with flange (solid, flexible or accordion), without built-in convexity, 4x4 inches or smaller, each
A4415	Ostomy skin barrier, with flange (solid, flexible or accordion), without built-in convexity, larger than 4x4 inches, each
A4416	Ostomy pouch, closed, with barrier attached, with filter (1-piece,) each
A4417	Ostomy pouch, closed, with barrier attached, with built-in convexity, with filter (1-piece), each
A4418	Ostomy pouch, closed; without barrier attached, with filter (1-piece) each
A4419	Ostomy pouch, closed; for use on barrier with non-locking flange, with filter (2-piece), each
A4421	Ostomy supply; miscellaneous
A4422	Ostomy pouch absorbent material (sheet, pad, crystal packet) for use in ostomy pouch to thicken liquid stomal output, each
A4423	Ostomy pouch, closed; for use on barrier with locking flange, with filter (2-piece), each
A4424	Ostomy pouch, drainable, with barrier attached, with filter (1-piece), each
A4425	Ostomy pouch, drainable; for use on barrier with non-locking flange, with filter (2-piece system), each
A4426	Ostomy pouch, drainable; for use on barrier with locking flange (2-piece system), each
A4427	Ostomy pouch, drainable; for use on barrier with locking flange, with filter (2 piece system), each
A4428	Ostomy pouch, urinary, with extended wear barrier attached, with faucet-type tap with valve (1-piece), each
A4429	Ostomy pouch, urinary, with barrier attached, with built-in convexity, with faucet-type tap with value (1-piece), each
A4430	Ostomy pouch, urinary, with extended wear barrier attached, with built-in convexity, with faucet-type tap with valve (1-piece), each
A4431	Ostomy pouch, urinary; with barrier attached, with faucet-type tap with valve (1- piece), each
A4432	Ostomy pouch, urinary; for use on barrier with non-locking flange, with faucet-type tap with valve (2-piece), each
A4433	Ostomy pouch, urinary; for use on barrier with locking flange (2-piece), each
A4434	Ostomy pouch, urinary; for use on barrier with locking flange, with faucet-type tap with valve (2-piece), each
A4435	Ostomy pouch, drainable, high output, with extended wear barrier (1-piece system), with or without filter, each
A4436	Irrigation supply; sleeve, reusable, per month
A4437	Irrigation supply; sleeve, disposable, per month

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services

HCPCS Code Reference	
A4455	Adhesive remover or solvent (for tape, cement or other adhesive), per ounce
A4456	Adhesive remover, wipes, any type, each
A5055	Stoma cap
A5056	Ostomy pouch, drainable, with extended wear barrier attached, with filter (1-piece), each
A5057	Ostomy pouch, drainable, with extended wear barrier attached, with built in convexity, with filter, (1-piece), each
A5061	Ostomy pouch, drainable; with barrier attached, (1-piece), each
A5062	Ostomy pouch, drainable; without barrier attached (1-piece), each
A5063	Ostomy pouch, drainable; for use on barrier with flange (2-piece system), each
A5071	Ostomy pouch, urinary; with barrier attached (1-piece), each
A5072	Ostomy pouch, urinary; without barrier attached (1-piece), each
A5073	Ostomy pouch, urinary; for use on barrier with flange (2-piece), each
A5093	Ostomy accessory; convex insert
A5120	Skin barrier, wipes or swabs, each
A5121	Skin barrier; solid, 6x6 or equivalent, each
A5122	Skin barrier; solid, 8x8 or equivalent, each

1. CMS. Local Coverage Determination (LCD) L33828 — Ostomy Supplies. Medicare Coverage Database. Centers for Medicare & Medicaid Services. Accessed January 2026.https://www.cms.gov/medicare-coverage-database/view/lcd.aspx?LCDId=33828

2. CMS. Local Coverage Determination (LCD) L33831 — Surgical Dressings. Medicare Coverage Database. Centers for Medicare & Medicaid Services. Accessed January 2026. https://www.cms.gov/medicare-coverage-database/view/lcd.aspx?LCDId=33831

Making life *easier*

Patient Education Resources
Scan here to enroll your patients into the Coloplast Care program

Healthcare Professional Login (coloplastcareconnect.us)

Scan here for ostomy patient education resources including How-To videos.

coloplastprofessional.us

Welcome to Coloplast
Supporting Products
Irrigation
2-Piece Click
2-Piece Flex
1-Piece
Fistula/Wound Management
Pediatric
Services